

Laboratory Refrigeration



Fan Circulated Laboratory Refrigeration
Advanced Laboratory Refrigeration
Sparkfree Laboratory Refrigeration

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Fan Circulated Laboratory Fridges

Fan circulated refrigerators are designed for temperature uniformity throughout the chamber of the refrigerator. Where not integral, all fan circulated fridges can be factory fitted with high/ low temperature alarms and digital temperature displays. Available in a range of sizes, all models feature a door lock for extra security and adjustable shelves so you can make the most of the space available.



[Labcold RLPF05042A Fan Air Circulated Laboratory Fridge 160 Litres](#)

Capacity: 160 Litre
Config: Under Bench
850mm H x 550mm W x 600mm D

Product Code: FPD-07541



[Labcold RLPF09043A Fan Air Circulated Laboratory Fridge 285 Litres](#)

Capacity: 285 Litre
Config: Free Standing
1435mm H x 600mm W x 625mm D

Product Code: FPD-07542



[Labcold RLHF13043 Fan Air Circulated Laboratory Fridge 390 Litres](#)

Capacity: 390 Litre
Config: Free Standing
1880mm H x 580mm W x 615mm D

Product Code: FPD-04922

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPF05042A Fan Air Circulated Laboratory Fridge 160 Litres



Product Code: FPD-07541

Dimensions: 850mm H x 550mm W x 600mm D

- Capacity 160 Litres
- Configuration: Underbench
- Fan assisted air circulation
- Door locks
- CFC, HCFC and Amonia free
- Low energy
- Auto defrosting
- Removable shelves
- Maximum Ambient Temperature up to 32°C
- 3 x Toughened Glass Shelves
- Operating Temperature +4°C
- Temperature Range: +2°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Alarm
- Temperature display

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPF09043A Fan Air Circulated Laboratory Fridge 285 Litres



Product Code: FPD-07542

Dimensions: 1435mm H x 600mm W x 625mm D

- Capacity 285 Litres
- Configuration: Free-standing
- Fan assisted air circulation
- Door locks
- CFC, HCFC and Amonia free
- Low energy
- Auto defrosting
- Removable shelves
- Maximum Ambient Temperature up to 32°C
- 5 Glass Shelves
- Operating Temperature +4°C
- Temperature Range: +2°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Alarm
- Temperature display

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLHF13043 Fan Air Circulated Laboratory Fridge 390 Litres



Product Code: FPD-04922

Dimensions: 1880mm H x 580mm W x 615mm D

- Capacity 390 Litres
- Configuration: Free-standing
- Fan assisted air circulation
- Door locks
- CFC, HCFC and Amonia free
- Temperature Display
- Auto defrosting
- Removable shelves
- Maximum Ambient Temperature up to 32°C
- 6 Shelves
- Operating Temperature +4°C
- Temperature Range: +2°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Alarm
- Temperature display

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Advanced Laboratory Fridges

Advanced fridges are completely constructed from powder coated mild steel with a stainless steel interior so are built to last. Insulated with Envirofoam® and fitted with micro processor controllers these fridges provide accurate temperature control with low GWP. All Advanced refrigerators have optimised air circulation for superior temperature control and feature a temperature display, high and low temperature alarms and door locks fitted as standard. For ease of positioning the larger cabinets have heavy duty wheels with a brake. The Advanced range of refrigerators can be fitted with access ports and remote alarm contacts if required



[Labcold RFR05043 Advanced Laboratory Fridge 150 Litres with Solid Door](#)

£1,905.00

Capacity: 150 Litre
Config: Under Bench
830mm H x 600mm W x 620mm D

Product Code: FPD-10127



[Labcold RFG05043 Advanced Laboratory Fridge 150 Litres with Glass Door](#)

£2,145.00

Capacity: 150 Litre
Config: Under Bench
830mm H x 600mm W x 620mm D

Product Code: FPD-10128



[Labcold RFR21042 Advanced Laboratory Fridge 577 Litres](#)

£3,650.00

Capacity: 577 Litre
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-04944



[Labcold RFG21042 Advanced Laboratory Fridge 577 Litres](#)

£4,210.00

Capacity: 577 Litre
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-04951



[Labcold RFR44042 Advanced Laboratory Fridge 1245 Litres](#)

£5,710.00

Capacity: 1245 Litre
Config: Free Standing
1980mm H x 1440mm W x 820mm D

Product Code: FPD-04945



[Labcold RFG44042 Advanced Laboratory Fridge 1245 Litres](#)

£6,685.00

Capacity: 1245 Litre
Config: Free Standing
1980mm H x 1440mm W x 820mm D

Product Code: FPD-04952

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR05043 Advanced Laboratory Fridge 150 Litres with Solid Door



Product Code: FPD-10127

Dimensions: 830mm H x 600mm W x 620mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Capacity: 150 litres (gross)
- Door Type: Solid
- Construction: AISI 304 grade stainless steel interior and exterior.
- Shelves Supplied: 2, Plastic coated wire mesh 530 x 325mm.
- Microprocessor temperature control with digital display to 1°C
- Hi/Lo temperature alarm - audible
- Fan circulated interior
- Door Lock
- CFC and HFC free
- Temperature set-point range: +2 to +8°C
- Automatic defrost with condensate removal
- Maximum ambient temperature: +43°C

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR05043 Advanced Laboratory Fridge 150 Litres with Solid Door



Product Code: FPD-10127

Dimensions: 830mm H x 600mm W x 620mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Capacity: 150 litres (gross)
- Door Type: Solid
- Construction: AISI 304 grade stainless steel interior and exterior.
- Shelves Supplied: 2, Plastic coated wire mesh 530 x 325mm.
- Microprocessor temperature control with digital display to 1°C
- Hi/Lo temperature alarm - audible
- Fan circulated interior
- Door Lock
- CFC and HFC free
- Temperature set-point range: +2 to +8°C
- Automatic defrost with condensate removal
- Maximum ambient temperature: +43°C

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR21042 Advanced Laboratory Fridge 577 Litres



Product Code: FPD-04944

Dimensions: 1980mm H x 720mm W x 820mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Configuration: Free-standing, solid door
- Capacity: 577 Litres
- Operating temperature: +4°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 4 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RAG21042 Advanced Laboratory Fridge 577 Litres



Product Code: FPD-04951

Dimensions: 1980mm H x 720mm W x 820mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Configuration: Free-standing, solid door
- Capacity: 577 Litres
- Operating temperature: +4°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 4 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RAFG44042 Advanced Laboratory Fridge 1245 Litres



Product Code: FPD-04952

Dimensions: 1980mm H x 1440mm W x 820mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Configuration: Free-standing, double glass door
- Capacity: 1245 Litres
- Operating temperature: +4°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 8 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RAFR44042 Advanced Laboratory Fridge 1245 Litres



Product Code: FPD-04945

Dimensions: 1980mm H x 1440mm W x 820mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Configuration: Free-standing, double solid door
- Capacity: 1245 Litres
- Operating temperature: +4°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 8 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Sparkfree Laboratory Fridges

The Sparkfree Fridges is engineered to store flammable materials safely with no sources of ignition inside the chambers. Everything that could conceivably cause a spark, such as the controller, is located outside of the refrigerator chamber making these fridges and freezers an important part of laboratory safety. Designed specifically for laboratories, these refrigerators are solidly built and come in a range of sizes, from a compact bench top to a large 400 litre refrigerator. All Sparkfree fridges are factory fitted with a door lock for optimum security.



[Labcold RLPR02042 Sparkfree Laboratory Fridge 60 Litres](#)

Capacity: 60 Litre
Config: Upright Refrigerator
572mm H x 485mm W x 520mm D
Product Code: FPD-04912



[LEC LR207C Laboratory Refrigerator 82 Litre](#)

Capacity: 82 Litre
Config: Bench Top
660mm H x 502mm W x 540mm D
Product Code: FPD-09794



[LEC LR307C Laboratory Refrigerator 107 Litre](#)

Capacity: 107 Litre
Config: Under Bench
830mm H x 495mm W x 620mm D
Product Code: FPD-09795



[Liebherr LKUexv 1610 Laboratory Refrigerator - 141 Litre](#)

Capacity: 141 Litre
Config: Bench Top
820mm H x 600mm W x 615mm D
Product Code: FPD-10260



[LEC LSR151 Laboratory Refrigerator 151 Litre](#)

Capacity: 151 Litre
Config: Under Bench
845mm H x 595mm W x 595mm D
Product Code: FPD-09796



[Labcold RLPR05042A Sparkfree Laboratory Fridge 160 Litres](#)

Capacity: 160 Litre
Config: Upright Refrigerator
850mm H x 550mm W x 600mm D
Product Code: FPD-07537



[LEC LSR288 Laboratory Refrigerator 288L Litre](#)

Capacity: 284 Litre
Config: Free Standing
1570mm H x 595mm W x 634mm D
Product Code: FPD-09797



[Labcold RLPR09043A Sparkfree Laboratory Fridge 285 Litres](#)

Capacity: 285 Litre
Config: Upright Refrigerator
1435mm H x 600mm W x 625mm D
Product Code: FPD-07538

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Sparkfree Laboratory Fridges

The Sparkfree Fridges is engineered to store flammable materials safely with no sources of ignition inside the chambers. Everything that could conceivably cause a spark, such as the controller, is located outside of the refrigerator chamber making these fridges and freezers an important part of laboratory safety. Designed specifically for laboratories, these refrigerators are solidly built and come in a range of sizes, from a compact bench top to a large 400 litre refrigerator. All Sparkfree fridges are factory fitted with a door lock for optimum security.



[Labcold RLHD13043 Sparkfree Laboratory Fridge 360 Litres](#)

Capacity: 360 Litre
Config: Upright Refrigerator
1750mm H x 585mm W x 655mm D
Product Code: FPD-04917



[Labcold RLSP1505 - 400 Litre Sparkfree Laboratory Fridge with Solid Door](#)

Capacity: 400 Litre
Config: Upright Refrigerator
1860mm H x 600mm W x 695mm D
Product Code: FPD-07707



[LEC LR1607C Laboratory Refrigerator 474 Litre](#)

Capacity: 474 Litre
Config: Free Standing
1688mm H x 695mm W x 910mm D
Product Code: FPD-09798

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPR02042 Sparkfree Laboratory Fridge 60 Litres



Product Code: FPD-04912

Dimensions: 572mm H x 485mm W x 520mm D

- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity 60 Litres
- Configuration: Counter top
- Sparkfree safety interiors
- CFC, HCFC and Amonia free
- Automatic defrosting
- Removable Shelves
- Maximum ambient temperatures up to 32°C
- Door Lock
- 2 Shelves
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temp Display
- Alarm
- Wall Mounting Kit

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LR207C Laboratory Refrigerator 82 Litre



Product Code: FPD-09794

Dimensions: 660mm H x 502mm W x 540mm D

Internal Dimensions: 540mm H x 415mm W x 340mm D

- +2°C to +10°C set temperature - Operates at +6°C for best possible storage of samples and chemicals
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- Small countertop size - For limited space needing biomedical and laboratory refrigeration
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid - Hinged, reversible
- Door Open Depth: 976mm
- Two wire Shelves: 2
- Shelf Size: 438x274x6mm
- Automatic Defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <43dB
- Insulation Thickness: 30mm
- Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 0.77

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LR307C Laboratory Refrigerator 107 Litre



Product Code: FPD-09795

Dimensions: 830mm H x 495mm W x 620mm D

Internal Dimensions: 540mm H x 403mm W x 400mm D

- +2°C to +10°C set temperature - Operates at +6°C for best possible storage of samples and chemicals
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark-free - There is nothing electrical inside the fridge so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Type: Solid - Hinged, reversible
- Door Open Depth: 1050mm
- Two Adjustable Front Feet
- Two Wire Shelves in Four Available Positions
- Shelf Size: 427x374x6mm
- Automatic Defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <43dB
- Insulation Thickness 35mm
- Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 0.86

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Liebherr LGUex 1500 Laboratory Freezer - 139 Litre



Product Code: FPD-10261

Dimensions: 820mm H x 600mm W x 615mm D

- Temperature range -9°C to -26°C
- Suitable for the storage of explosive or highly inflammable substances
- Digital Temperature Control & Display
- Temperature deviation alarm with volt free contacts
- Door open alarm
- Manual Defrost
- Door Lock and 2 keys
- CFC & Ammonia free
- 3 x drawers + 1 Basket
- 220-230v 50Hz
- Max. Ambient +38°C

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Liebherr LKUexv 1610 Laboratory Refrigerator - 141 Litre



Product Code: FPD-10260

Dimensions: 820mm H x 600mm W x 615mm D

- Temperature range +3°C to +8°C
- Suitable for the storage of explosive or highly inflammable substances
- Digital Temperature Control & Display
- Temperature deviation alarm with volt free contacts
- Door open alarm
- Auto Defrost
- Door Lock and 2 keys
- CFC & Ammonia free
- 3 x shelves
- 220-230v 50Hz
- Max. Ambient +38°C
- 2 Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

LEC LSR151 Laboratory Refrigerator 151 Litre



Product Code: FPD-09796

Dimensions: 845mm H x 595mm W x 595mm D

Internal Dimensions: 705mm H x 500mm W x 450mm D

- Set temperature - Operates at +X°C for best possible storage of samples and chemicals
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: +2°C to +10°C
- Colour: white/green
- Door Type: Solid - Hinged, reversible
- Open Door Sensor: Yes
- Lock
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Temperature Display on Panel
- Temperature Alarm
- Three Wire Shelves Six Positions
- Shelf Size: 320x520x6mm
- Two Adjustable Door Wire Bins
- Automatic Defrost
- External Digital Controller
- Climate Class: SN
- Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Insulation Thickness: 45mm
- Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPR05042A Sparkfree Laboratory Fridge 160 Litres



Product Code: FPD-07537

Dimensions: 850mm H x 550mm W x 600mm D

- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity 160 Litres
- Configuration: Underbench lockable
- Sparkfree safety interiors
- CFC, HCFC and Amonia free
- Automatic defrosting
- Removable Shelves
- Maximum ambient temperatures up to 32°C
- Door Lock
- Low Energy
- 4 Toughened Glass Shelves
- 1 x Drawer
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temp Display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSR288 Laboratory Refrigerator 288L Litre



Product Code: FPD-09797

Dimensions: 1570mm H x 595mm W x 634mm D

Internal Dimensions: 1382mm H x 480mm W x 440mm D

- Set temperature - Operates at +X°C for best possible storage of samples and chemicals
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- External Finish: white
- Internal Finish: white
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: +2°C to +10°C
- Colour: White/green
- Door Type: Solid - Hinged, reversible
- Door Open Depth: 1185mm
- Open Door Sensor
- Lock
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Temperature Display on Panel
- Temperature Alarm
- Five Wire Shelves, 12 Positions
- Shelf Size: 540x497x6mm
- One Solid Moulded Drawer
- Four Solid Adjustable Door Bins
- Automatic Defrost
- External Digital Controller
- Climate Class: SN
- Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Insulation Thickness: 60mm
- Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPR09043A Sparkfree Laboratory Fridge 285 Litres



Product Code: FPD-07538

Dimensions: 1435mm H x 600mm W x 625mm D

- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity 285 Litres
- Configuration: Free-standing lockable
- Sparkfree safety interiors
- CFC, HCFC and Amonia free
- Automatic defrosting
- Removable Shelves
- Maximum ambient temperatures up to 32°C
- Low Energy
- Door Lock
- 5 Glass Shelves
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temp Display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLHD13043 Sparkfree Laboratory Fridge 360 Litres



Product Code: FPD-04917

Dimensions: 1750mm H x 585mm W x 655mm D

- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity 360 Litres
- Configuration: Free-standing
- Sparkfree safety interiors
- CFC, HCFC and Amonia free
- Automatic defrosting
- Removable Shelves
- Maximum ambient temperatures up to 32°C
- 6 Shelves
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temp Display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLSP1505 - 400 Litre Sparkfree Laboratory Fridge with Solid Door



Product Code: FPD-07707

Dimensions: 1860mm H x 600mm W x 695mm D

- A purpose-built refrigerator for laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future. Fan circulation of chamber air provides much improved temperature uniformity.
- Capacity: 400 litres (gross)
- Door Type: Solid
- Construction: White coated steel interior and exterior
- Shelves Supplied: 6, Plastic coated wire mesh
- Sparkfree Interior
- Fan assisted air circulation powered by external motor
- Digital display of temperature
- Audible alarm with mute facility
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Door Lock
- Easy-grip door handle
- Access Port
- CFC, HCFC and Ammonia free
- Accurate temperature
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LR1607C Laboratory Refrigerator 474 Litre



Product Code: FPD-09798

Dimensions: 1688mm H x 695mm W x 910mm D

Internal Dimensions: 1170mm H x 550mm W x 670mm D

- +2°C to +10°C set temperature - Operates at +6°C for best possible storage of samples and chemicals
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Front lockable roller castors - To make this large fridge easier to move for cleaning and maintenance
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Type: Solid - Hinged, R/H only
- Door Open Depth: 1596mm
- Lock
- Adjustable Feet: Castors
- Five Wire Shelves 30 Positions: 30
- Shelf Size: 536x609x6mm
- Automatic Defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Insulation Thickness: 65mm
- Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Delivery

48 Hour delivery is available on most refrigerators through-out the UK mainland, larger units may require a longer lead time, please contact our sales staff whilst ordering.

Approximate Delivery lead times can be confirmed during our quotation process, typical delivery time is 10-14 days from receipt of payment.

**Drivers will not carry goods into buildings.

Goods normally arrive by truck on a pallet and are delivered Kerbside. Please note you may need lifting equipment or a number of strong workers to handle the order when it arrives. If you have any special requirements please [contact us](#).

Report any damages or shortages immediately, please read delivery notes carefully and use contact numbers provided or call FPD Group Ltd.

For more detail refer to [FPD Group Ltd Terms and Conditions](#) documentation.

**Large Products will require a two person delivery and installation, this is built into the price, the larger fridges also may require a longer delivery time of 3-4 weeks please check whilst ordering.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP[™]

Pharmacy Fridge Exports

Since 1961 FPD Group Ltd's family have been supplying medical equipment. Our extensive experience and technical knowledge now concentrates on providing clients with a quality selection of products for hospitals and healthcare organisations worldwide. We receive export enquiries and orders from all parts of the world and have established customers in UK, Europe, Africa and Asia. Export business activities are supervised by either our UK or Asia offices.

Our client base ranges from overseas ministry of health, defence, hospitals, surgeries, clinics and doctors to aid agencies and disaster relief. FPD Group Ltd export almost any medical product listed on our websites.

FPD estimators have the expertise to supply a single item or to fully equip many types of healthcare facilities. As well as the initial supply of equipment we also offer room design services using our in-house Autodesk Revit Architecture 2013 CAD systems.

** Export items are normally subject to minimum order values, additional protective packing charges and export shipping costs. Prices shown are Ex Works unless otherwise agreed in writing.*



FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

About Pharmacy Fridge

Pharmacy Fridge supply a wide range of fridges meeting the specifications of the healthcare industry. To buy simply add your required product to the shopping cart and submit your order.

For dispensary, surgery, ward and general use for the storage of drugs and vaccines under controlled conditions in line with the "Code of Ethics and Standards" of the Royal Pharmaceutical Society of Great Britain. Where necessary to achieve the required temperature control internal air-circulating fans are fitted.

Temperature display, alarms, max/min temperature recording, and door locks are standard and all models are CFC, HCFC and Ammonia free.

Are you storing your vaccines correctly? Recent research has shown an alarming level of use of domestic refrigerators to store vaccines. Did you know that not only can this damage the effectiveness of the vaccines it also poses a security risk for medicines that should be locked away?

FPD Group Ltd

FPD Group Limited is a family business with over 50 years experience in the supply and installation of a wide range of Healthcare Equipment, Medical Refrigeration and Hygienic Storage Furniture. Our customer base includes organisations such as NHS Trusts, PCT's, HMP Prison service, Ministry of Defence, GP Practices, Pharmacy Shops and Dispensaries, Universities, Schools, Dental Practices, Laboratories, Veterinary Practices, Care Homes and Private Buyers. We also export to the Irish Republic, Asia, United Arab Emirates and Qatar. In the UK the Company is based in Exeter, Devon with an additional office in Asia.

Our product range includes pharmacy storage solutions, refrigeration, controlled drugs storage, medical trolleys, clinical waste bins, medical lighting solutions, medical and healthcare furniture ranges and a number of bespoke solutions on request.

FPD Group has an in-house CAD visualisation department and accepts most file formats from your architects' plans, please ask for a quotation.

We are fully committed to the provision of a comprehensive range of professional products, services and furniture. Our philosophy is to provide a complete and appropriate solution to our customers' requirements.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Contact Pharmacy Fridge

FPD Group Ltd have over 40 years experience in the healthcare and shopfitting industry. We offer consultancy, design and installation services. Along with our Pharmacy Shelving and Storage solutions we also provide a wide range of medical equipment and healthcare furniture. Please visit www.fpd-group.co.uk for more information on our medical products.

For general enquiries complete the form below and we will contact you ASAP, please provide as much information as possible including dimensions where required. Please email any plans or photos of installations directly to sales@medicine-fridge.co.uk.

| | | | |
|-----------|----------------------|------------|----------------------|
| Name | <input type="text"/> | Address | <input type="text"/> |
| Company | <input type="text"/> | ... | <input type="text"/> |
| Telephone | <input type="text"/> | City | <input type="text"/> |
| Fax | <input type="text"/> | State | <input type="text"/> |
| Email | <input type="text"/> | PostalCode | <input type="text"/> |
| | | Country | <input type="text"/> |
| Comments | <input type="text"/> | | |

Submit Enquiry

Sales

Tel: 0113 350 8696
Fax: 0115 727 0123
Email: sales@pharmacy-shelving.co.uk

Payments & Billing

Tel: 01392 927 567

Address

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

How To Order

NHS Purchase Order, MOD Purchase Order, Local Authority Purchase Order

All official government, NHS, MOD, Police, Educational and Local Authority purchase orders are accepted. Fax your purchase order to 0115 727 0123 or email to info@pharmacy-fridge.co.uk.

Bank Transfer

Sent bank transfers to:

| | | | |
|----------------|-------------------|-------|------------------------|
| Payee | FPD Group Limited | | |
| Bank | HSBC Bank plc | | |
| Sortcode | 401903 | IBAN | GB97MIDL40190302074613 |
| Account Number | 02074613 | Swift | MIDLGB2104A |

If possible send confirmation of your transaction, to speed up the order process.

Cheques

Make cheques payable to FPD Group Limited and post to:

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Credit/Debit Card

Payment from UK based businesses can be taken over the telephone or processed via the website [shopping cart](#). Customers with issued quotations may process their order online using our [Purchase Code](#) web page.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk