

Medical Refrigeration



Pharmacy Refrigeration
Advanced Pharmacy Refrigeration
Laboratory Refrigeration
Low Temperature Refrigeration
Portable Refrigeration
Vaccine Carriers
Ward Refrigeration
Bloodbanks
Laboratory Incubators
Icemakers & Water Dispensers

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Pharmacy & Vaccine Refrigeration



Pharmacy Refrigeration
LEC Medical
Labcold
Shoreline Medical
Dometic Silenco
Advanced Pharmacy Refrigeration

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[28 Litre Silenco Lockable Pharmacy Absorption Fridge](#)

£430.00

Capacity: 28 Litre
Config: Bench Top
580mm H x 422mm W x 393mm D

Product Code: FPD-09818



[Labcold IntelliCold™ RLDF0110 - 36 Litre Pharmacy Fridge with Solid Door](#)

£495.00

Capacity: 36 Litre
Config: Bench Top/Wall Mounted
535mm H x 445mm W x 465mm D

Product Code: FPD-09835



[LEC PG102C - Pharmacy Fridge 41 Litre with Glass Door](#)

£695.00

Capacity: 41 Litre
Config: Bench Top
702mm H x 379mm W x 500mm D

Product Code: FPD-09752



[LEC PE102C - Pharmacy Fridge 41 Litre with Solid Door](#)

£585.00

Capacity: 41 Litre
Config: Bench Top
702mm H x 379mm W x 500mm D

Product Code: FPD-09742



[LEC PE109C - Pharmacy Fridge 45 Litre with Solid Door](#)

£441.00

Capacity: 45 Litre
Config: Bench Top
530mm H x 470mm W x 465mm D

Product Code: FPD-09743



[Shoreline SM45 - 45 Litres Pharmacy Fridge with Solid Door](#)

£478.00

Capacity: 45 Litre
Config: Bench Top
715mm H x 380mm W x 375mm D

Product Code: FPD-06627



[53 Litre Silenco Lockable Pharmacy Absorption Fridge](#)

£515.00

Capacity: 53 Litre
Config: Bench Top
592mm H x 486mm W x 494mm D

Product Code: FPD-09819



[Shoreline SM60G - 60 Litres Pharmacy Fridge with Glass Door](#)

£718.00

Capacity: 60 Litre
Config: Bench Top
660mm H x 470mm W x 455mm D

Product Code: FPD-04358

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[Labcold IntelliCold™ RLDF0210 - 66 Litre Pharmacy Fridge with Solid Door](#)

£590.00

Capacity: 66 Litre
Config: Bench Top/Wall Mounted
725mm H x 450mm W x 485mm D
Product Code: FPD-09836



[Shoreline SM75 - 75 Litres Pharmacy Fridge with Solid Door](#)

£608.00

Capacity: 75 Litre
Config: Bench Top
790mm H x 450mm W x 430mm D
Product Code: FPD-07731



[LEC PE207C - Pharmacy Fridge 82 Litre with Solid Door](#)

£754.00

Capacity: 82 Litre
Config: Bench Top
660mm H x 502mm W x 540mm D
Product Code: FPD-09744



[LEC PG207C - Pharmacy Fridge 82 Litre with Glass Door](#)

£999.00

Capacity: 82 Litre
Config: Bench Top
660mm H x 502mm W x 450mm D
Product Code: FPD-09753



[LEC PG307C - Pharmacy Fridge 107 Litre with Glass Door](#)

£1,213.00

Capacity: 107 Litre
Config: Under Bench
830mm H x 495mm W x 620mm D
Product Code: FPD-09754



[LEC PE307C - Pharmacy Fridge 107 Litre with Solid Door](#)

£980.00

Capacity: 107 Litre
Config: Under Bench
830mm H x 495mm W x 620mm D
Product Code: FPD-09745



[Shoreline SM161G - 141 Litres Pharmacy Fridge with Glass Door](#)

£1,128.00

Capacity: 141 Litre
Config: Under Bench
835mm H x 600mm W x 615mm D
Product Code: FPD-08130



[Shoreline SM161 - 141 Litres Pharmacy Fridge with Solid Door](#)

£938.00

Capacity: 141 Litre
Config: Under Bench
835mm H x 600mm W x 615mm D
Product Code: FPD-08131

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[Shoreline SM150 - 150 Litres Pharmacy Fridge with Solid Door](#)

£798.00

Capacity: 150 Litre
Config: Under Bench
855mm H x 550mm W x 600mm D

Product Code: FPD-04361



[Labcold IntelliCold™ RLDF0510 - 150 Litre Pharmacy Fridge with Solid Door](#)

£860.00

Capacity: 150 Litre
Config: Under Bench
820mm H x 600mm W x 610mm D

Product Code: FPD-09837



[Labcold IntelliCold™ RLDG0510 - 150 Litre Pharmacy Fridge with Glass Door](#)

£1,060.00

Capacity: 150 Litre
Config: Under Bench
820mm H x 600mm W x 610mm D

Product Code: FPD-09838



[LEC PSR151UK - Pharmacy Fridge 151 Litre with Solid Door](#)

£895.00

Capacity: 151 Litre
Config: Under Bench
845mm H x 595mm W x 635mm D

Product Code: FPD-09746



[LEC PGR151UK - Pharmacy Fridge 151 Litre with Glass Door](#)

£1,195.00

Capacity: 151 Litre
Config: Under Bench
845mm H x 595mm W x 635mm D

Product Code: FPD-09755



[Shoreline SM260 - 260 Litres Pharmacy Fridge with Solid Door](#)

£1,098.00

Capacity: 260 Litre
Config: Free Standing
1225mm H x 600mm W x 600mm D

Product Code: FPD-04363



[Shoreline SM260G - 260 Litres Pharmacy Fridge with Glass Door](#)

£1,348.00

Capacity: 260 Litre
Config: Free Standing
1225mm H x 600mm W x 600mm D

Product Code: FPD-04359

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[LEC PSR273UK - Pharmacy Fridge 273
Litre with Solid Door](#)

£1,195.00

Capacity: 273 Litre
Config: Free Standing
1565mm H x 595mm W x 670mm D
Product Code: FPD-09747



[LEC PGR273UK - Pharmacy Fridge 273
Litre with Glass Door](#)

£1,595.00

Capacity: 273 Litre
Config: Free Standing
1565mm H x 595mm W x 660mm D
Product Code: FPD-09756



[Labcold IntelliCold™ RLDG1010 - 300
Litre Pharmacy Fridge with Solid Door](#)

£1,160.00

Capacity: 300 Litre
Config: Free Standing
1500mm H x 600mm W x 660mm D
Product Code: FPD-09839



[Labcold IntelliCold™ RLDG1010 - 300
Litre Pharmacy Fridge with Glass Door](#)

£1,510.00

Capacity: 300 Litre
Config: Free Standing
1500mm H x 600mm W x 600mm D
Product Code: FPD-09840



[LEC PGR334UK - Pharmacy Fridge 334
Litre with Glass Door](#)

£1,995.00

Capacity: 334 Litre
Config: Free Standing
1860mm H x 595mm W x 660mm D
Product Code: FPD-09757



[LEC PSR334UK - Pharmacy Fridge 334
Litre with Solid Door](#)

£1,795.00

Capacity: 334 Litre
Config: Free Standing
1860mm H x 595mm W x 670mm D
Product Code: FPD-09748



[Shoreline SM360G - 360 Litres
Pharmacy Fridge with Glass Door](#)

£1,668.00

Capacity: 360 Litre
Config: Free Standing
1595mm H x 600mm W x 600mm D
Product Code: FPD-04356



[Shoreline SM360 - 360 Litres Pharmacy
Fridge with Solid Door](#)

£1,398.00

Capacity: 360 Litre
Config: Free Standing
1595mm H x 600mm W x 600mm D
Product Code: FPD-04364

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[Labcold IntelliCold™ RLDF1510 - 430 Litre Pharmacy Fridge with Solid Door](#)

£1,960.00

Capacity: 430 Litre
Config: Free Standing
1860mm H x 600mm W x 700mm D

Product Code: FPD-09843



[Labcold IntelliCold™ RLDG1510 - 430 Litre Pharmacy Fridge with Glass Door](#)

£2,125.00

Capacity: 430 Litre
Config: Free Standing
1860mm H x 600mm W x 700mm D

Product Code: FPD-09844



[LEC PE1607C - Pharmacy Fridge 444 Litre with Solid Door](#)

£2,975.00

Capacity: 444 Litre
Config: Free Standing
1688mm H x 685mm W x 910mm D

Product Code: FPD-09749



[LEC PG1607C - Pharmacy Fridge 444 Litre with Glass Door](#)

£3,229.00

Capacity: 444 Litre
Config: Free Standing
1688mm H x 685mm W x 910mm D

Product Code: FPD-09758



[Shoreline SM541 - 544 Litres Pharmacy Fridge with Solid Door](#)

£2,485.00

Capacity: 544 Litre
Config: Free Standing
1745mm H x 750mm W x 700mm D

Product Code: FPD-07642



[Shoreline SM541G - 544 Litres Pharmacy Fridge with Glass Door](#)

£2,850.00

Capacity: 544 Litre
Config: Free Standing
1745mm H x 750mm W x 700mm D

Product Code: FPD-07643



[LEC PSR600UK - Pharmacy Fridge 600 Litre with Solid Door](#)

£3,500.00

Capacity: 600 Litre
Config: Free Standing
1995mm H x 690mm W x 840mm D

Product Code: FPD-10199



[LEC PGR600UK - Pharmacy Fridge 600 Litre with Glass Door](#)

£4,000.00

Capacity: 600 Litre
Config: Free Standing
1995mm H x 690mm W x 840mm D

Product Code: FPD-10201

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Pharmacy Fridges

All pharmacy fridges are developed with reliability and ease of use in mind. Our pharmacy refrigeration range conforms to the guidelines set out by the RPSGB and maintains an internal temperature of between +2°C and +8°C. We have added a number of improvements to certain popular models to ensure we are offering you the best possible secure and reliable refrigerated storage. Many of our fridges are available with different plug options for Europe and Australasia, or with compressors operating at 60Hz.



[LEC PSR1200UK - Pharmacy Fridge](#)
[1200 Litre with Solid Door](#)

£6,000.00

Capacity: 1200 Litre
Config: Free Standing
1995mm H x 1330mm W x 840mm D
Product Code: FPD-10200



[LEC PGR1200UK - Pharmacy Fridge](#)
[1200 Litre with Glass Door](#)

£6,500.00

Capacity: 1200 Litre
Config: Free Standing
1995mm H x 1330mm W x 840mm D
Product Code: FPD-10202

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

28 Litre Silenco Lockable Pharmacy Absorption Fridge



Product Code: FPD-09818

Dimensions: 580mm H x 422mm W x 393mm D

- The perfect solution for hospitals and health centres to keep heat sensitive medicines between 2-8°C. With magnetic sealing door closing. LED light and acoustic signal as standard. Signal indicates that the door has been opened too long. Signal starts 20 seconds after door opening.
- Internal fan to equalize cabinet temperature and avoid frost formation on the evaporators fin. Trays with transparent front board to keep the inside from warming up by blocking hot air penetration into the cabinet when opening the door. It also displays the min and max temperatures from the last reset.
- Digital thermometer to display the actual inside temperature
- Outside jacket material: polystyrol
- Cooling unit control: thermofuse
- Colour: white
-

Delivery

Carriage

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold IntelliCold™ RLDF0110 - 36 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09835

Dimensions: 535mm H x 445mm W x 465mm D

- Powered by the unique IntelliCold™, this compact pharmacy and vaccine fridge has been designed from the ground up for the safe storage of cold chain medicines.
- Featuring an all metal interior and exterior, it is built to meet the demands of a professional setting and comes with all the features you would expect on a much larger and expensive professional pharmacy refrigerator.
- Small enough to sit on a bench top or to be wall mounted by way of integral brackets, this versatile pharmacy fridge is ideal for surgery situations, care homes or schools, indeed anywhere where there is a requirement to store limited amount of temperature sensitive vaccines and pharmaceuticals in optimum conditions.
- Configuration: Wall mounted/Bench top
- Operating temp (nominal) 5°C
- Temperature range 2 to 8°C
- Door Type: Solid
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up
- Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Two Shelves

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PG102C - Pharmacy Fridge 41 Litre with Glass Door



Product Code: FPD-09752

Dimensions: 702mm H x 379mm W x 500mm D

Internal Dimensions: 400mm H x 245mm W x 310mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Small countertop size- For limited space needing medical refrigeration
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-R/H only, Door Open Depth: 810mm
- Two wire shelves, Number of Shelf Positions: 5, Size: 270x308x6mm
- Automatic defrost
- Climate Class: SN
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <45dB
- Insulation Thickness: 40mm, Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.68

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PE102C - Pharmacy Fridge 41 Litre with Solid Door



Product Code: FPD-09742

Dimensions: 702mm H x 379mm W x 500mm D

Internal Dimensions: 400mm H x 245mm W x 310mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Small countertop size - For limited space needing medical refrigeration
- Comes supplied with Wall Bracket
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Door Finish: Painted steel
- Colour: White/green
- Door Type: Solid, Hinged R/H, Door Open Depth: 810Rmm
- Two wire shelves, available positions, Shelf Size: 270x308x6mm
- Automatic defrost
- Climate Class: SN
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <45dB
- Insulation Thickness: 40mm, Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.56

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PE109C - Pharmacy Fridge 45 Litre with Solid Door



Product Code: FPD-09743

Dimensions: 530mm H x 470mm W x 465mm D

Internal Dimensions: 400mm H x 410mm W x 270mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- Small countertop size - For limited space needing medical refrigeration
- Comes supplied with Wall Bracket
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 910mm
- Two Adjustable Feet at front
- One wire shelf, Size: 402x200x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <45dB
- Insulation Thickness: 30mm, Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.68

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM45 - 45 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-06627

Dimensions: 715mm H x 380mm W x 375mm D

Internal Dimensions: 445mm H x 300mm W x 250mm D

- 45 litres Capacity
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and 2 keys
- CFC & Ammonia free
- Adjustable Shelving x 2
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- 220-230v 50Hz
- Max. Ambient +32°C
- Prices Includes Delivery Mainland UK
- 2 Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP™

53 Litre Silenco Lockable Pharmacy Absorption Fridge



Product Code: FPD-09819

Dimensions: 592mm H x 486mm W x 494mm D

- Digital thermometer to display the actual inside temperature. It also displays the min and max temperatures from the last reset. LED light and acoustic signal as standard. Signal indicates that the door has been opened too long. Signal starts 20 seconds after door opening.
- Internal fan to equalize cabinet temperature and avoid frost formation on the evaporators fin. Trays with transparent front board to keep the inside from warming up by blocking hot air penetration into the cabinet when opening the door.
- Lock/door closing: magnetic sealing
- Cabinet temperature: 2°C to 8°C
- Outside jacket material: polystyrol
- Cooling unit control: thermofuse
- Colour: white

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM60G - 60 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-04358

Dimensions: 660mm H x 470mm W x 455mm D

Internal Dimensions: 525mm H x 375mm W x 310mm D

- Glass Door
- Capacity 60 Litres
- 2 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Labcold IntelliCold™ RLDF0210 - 66 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09836

Dimensions: 725mm H x 450mm W x 485mm D

- Powered by the unique IntelliCold™, this compact pharmacy and vaccine fridge has been designed from the ground up for the safe storage of cold chain medicines.
- Featuring an all metal interior and exterior, it is built to meet the demands of a professional setting and comes with all the features you would expect on a much larger and expensive professional pharmacy refrigerator.
- Small enough to sit on a bench top or to be wall mounted by way of integral brackets, this versatile pharmacy fridge is ideal for surgery situations, care homes or schools. Fully lockable and with the added extra of automatic 24/7 integral temperature logging, this pharmacy fridge is designed to be compliant with all current professional guidelines including the DH Green Book, NaTHNaC, GPhC/RPS and RCVS/BSVA.
- Configuration: Bench top/wall mounted
- Door Type: Solid
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8deg;C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Three Shelves



Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM75 - 75 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-07731

Dimensions: 790mm H x 450mm W x 430mm D

Internal Dimensions: 525mm H x 375mm W x 305mm D

- 75 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and 2 keys
- CFC & Ammonia free
- Adjustable Shelving x 2
- Meets RPSGB, NHS, WHO, RCVS performance criteria
- 220-240v 50/60Hz
- Max. Ambient +32°C
- Made in England
- 2 Years Part & Labour Warranty (UK only)
- Price Includes Delivery to Mainland UK

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

LEC PE207C - Pharmacy Fridge 82 Litre with Solid Door



Product Code: FPD-09744

Dimensions: 660mm H x 502mm W x 540mm D

Internal Dimensions: 540mm H x 419mm W x 315mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Door Finish: Painted steel
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 976mm
- Two wire shelves, Size: 438x274x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <43dB
- Insulation Thickness: 30mm
- Energy Consumption kw/24h: 1.46

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PG207C - Pharmacy Fridge 82 Litre with Glass Door



Product Code: FPD-09753

Dimensions: 660mm H x 502mm W x 450mm D

Internal Dimensions: 540mm H x 419mm W x 320mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-reversible, Door Open Depth: 976mm, Replacement Door Seal
- Open Door Sensor
- Two wire Shelves, Number of Shelf Positions: 2, Size: 438x274x6mm
- Interior Light (Medical)
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level: <43dB
- Insulation Thickness: 30mm
- Energy Consumption kw/24h: 1.58

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PG307C - Pharmacy Fridge 107 Litre with Glass Door



Product Code: FPD-09754

Dimensions: 830mm H x 495mm W x 620mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-reversible, Door Open Depth: 1050mm, Replacement Door Seal
- Adjustable Feet: 1, front
- Two wire shelves, Number of Shelf Positions: 4, Size - 427x374x6mm
- Interior Light (Medical)
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230, Aux Output
- Insulation Thickness: 35mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PE307C - Pharmacy Fridge 107 Litre with Solid Door



Product Code: FPD-09745

Dimensions: 830mm H x 495mm W x 620mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted steel
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1050mm
- Adjustable Feet: 1, front
- Two wire shelves, Number of Shelf Positions: 4, size: 427x374x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 35mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM161 - 141 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-08131

Dimensions: 835mm H x 600mm W x 615mm D

Internal Dimensions: 667mm H x 441mm W x 439mm D

- Solid Door
- 141 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display with ShoreTouch? controls
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Self-Closing Door with Lock and 2 keys
- High Efficiency - Thick Wall Insulation
- Adjustable Shelving x 3
- Interior Light
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +43°C
- Made in England
- 2 Years Parts & Labour Warranty (UK only)
- Supplied with Certificate of Conformity

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM161G - 141 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-08130

Dimensions: 835mm H x 600mm W x 615mm D

Internal Dimensions: 667mm H x 441mm W x 439mm D

- Glass Door
- 141 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display with ShoreTouch controls
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Self-Closing Door with Lock and 2 keys
- High Efficiency - Thick Wall Insulation
- Adjustable Shelving x 3
- Interior Light
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +43°C
- Made in England
- 2 Years Parts & Labour Warranty (UK only)
- Supplied with Certificate of Conformity

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Shoreline SM260 - 260 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-04363

Dimensions: 1225mm H x 600mm W x 600mm D

Internal Dimensions: 1047mm H x 513mm W x 441mm D

- Capacity 260 Litres Gross, 247 Litres Net
- 4 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48
HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP_{TM}

Shoreline SM150 - 150 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-04361

Dimensions: 855mm H x 550mm W x 600mm D

Internal Dimensions: 665mm H x 460mm W x 385mm D

- Capacity 150 Litres
- 3 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP_{TM}

Labcold IntelliCold™ RLDG0510 - 150 Litre Pharmacy Fridge with Glass Door



Product Code: FPD-09838

Dimensions: 820mm H x 600mm W x 610mm D



- A professional 150 litre under bench pharmacy refrigerator powered by the Labcold IntelliCold™ controller. Designed for the safe storage of cold chain pharmaceuticals, this refrigerator assists busy professionals meet the regulatory requirements with an easy to operate intuitive four button controller and 13 month temperature data storage which can be downloaded to any PC or Mac. For ease of viewing the contents and to minimize door openings, the light can be set to operate independently of the door and even left on all the time.
- Reliable and robust, this refrigerator is engineered to meet the demands of a busy practice, surgery or pharmacy and is backed by the Labcold FREE 5 year warranty.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Underbench
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Three Shelves
- Door Type: Glass
- Access port

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PSR273UK - Pharmacy Fridge 273 Litre with Solid Door



Product Code: FPD-09747

Dimensions: 1565mm H x 595mm W x 670mm D

Internal Dimensions: 1382mm H x 480mm W x 445mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Anti-microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Interior Xw light - For better visibility of contents
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1185mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Five wire Shelves, Number of Shelf Positions: 12
- Interior Light (Medical)
- Automatic defrost
- Climate Class: SN
- Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 220-240v/50hz, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold IntelliCold™ RLDF0510 - 150 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09837

Dimensions: 820mm H x 600mm W x 610mm D



- A professional 150 litre under bench pharmacy refrigerator powered by the Labcold IntelliCold™ controller. Designed for the safe storage of cold chain pharmaceuticals, this refrigerator helps busy professionals meet regulatory requirements laid down in the DH Green Book with an easy to operate intuitive four button controller and 13 month temperature data storage. This data can be downloaded to a free SD card (available when the warranty is registered) and stored on a PC or Mac. Reliable and robust, this refrigerator is engineered to meet the demands of a busy practice or pharmacy and is backed by the Labcold FREE 5 year parts and labour warranty.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Underbench
- Door Type: Solid
- Operating temp (nominal): 5°C
- Temperature range: 2 - 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Three Shelves
- Access port

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR151UK - Pharmacy Fridge 151 Litre with Glass Door



Product Code: FPD-09755

Dimensions: 845mm H x 595mm W x 635mm D

Internal Dimensions: 705mm H x 500mm W x 450mm D

- PSR151 solid door refrigerator is ideal for the under-counter storage of vaccines and comes complete with our new advanced controller and a fitted lock for added security
- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Antimicrobial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening.
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted steel
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1165mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three wire shelves, Number of Shelf Positions: 6, Size: 400x520x6mm
- Automatic defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 45mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR273UK - Pharmacy Fridge 273 Litre with Glass Door



Product Code: FPD-09756

Dimensions: 1565mm H x 595mm W x 660mm D

Internal Dimensions: 1382mm H x 480mm W x 445mm D

- PGR273 glass door refrigerator is ideal for the storage of vaccines and comes complete with our new advanced controller and a fitted lock for added security.
- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Interior Xw light - For better visibility of contents without having to open the door
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted aluminium
- Door Type: Glass, Hinged-reversible, Door Open Depth: 1180mm, Replacement Door Seal
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Five wire Shelves, Number of Shelf Positions: 12
- Automatic defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 220-240v/50hz, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PSR151UK - Pharmacy Fridge 151 Litre with Solid Door



Product Code: FPD-09746

Dimensions: 845mm H x 595mm W x 635mm D

Internal Dimensions: 705mm H x 500mm W x 450mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Antimicrobial handle
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening.
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted steel
- Colour: White/green
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1165mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three wire Shelves, Number of Shelf Positions: 6, size: 400x520x6mm
- Automatic defrost
- Climate Class: SN
- Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 45mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM260G - 260 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-04359

Dimensions: 1225mm H x 600mm W x 600mm D

Internal Dimensions: 1047mm H x 513mm W x 441mm D

- Glass Door
- Interior Light
- Capacity 260 Litres Gross, 252 Litres Net
- 4 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold IntelliCold™ RLDG1010 - 300 Litre Pharmacy Fridge with Glass Door



Product Code: FPD-09840

Dimensions: 1500mm H x 600mm W x 600mm D

- With a toughened safety glass door, this 300 litre pharmacy fridge, fitted with the Labcold IntelliCold™ controller, is purposed designed for the safe storage of temperature sensitive vaccines and medicines. With an external light switch that can operate independently of the fridge door, this model can not only form an attractive display feature but allows for quick visual stock takes without the need for opening the door. Backed by the Labcold FREE 5 year parts and labour warranty, there has never been a better way to store cold chain pharmaceuticals.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Upright
- Door Type: Glass
- Operating temp (nominal): 5°C
- Temperature range: 2 - 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Six Shelves
- Access port

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.



FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold IntelliCold™ RLDF1010 - 300 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09839

Dimensions: 1500mm H x 600mm W x 660mm D

- A freestanding 300 litre pharmacy refrigerator fitted with the Labcold IntelliCold™ controller to ensure the safe storage of cold chain pharmaceuticals. The easy to operate push button controller clearly displays the internal temperature and displays minimum and maximum temperature since the last temperature re-set at the touch of the button. Immediate high/low temperature alarms and controller battery back up which means it functions even in the event of a power failure, makes this fridge one of the safest for vaccine storage. For extra security, this refrigerator is backed by the Labcold FREE 5 year parts and labour warranty.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Upright
- Door Type: Solid
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up Door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Six Shelves
- Access port



Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PG1607C - Pharmacy Fridge 444 Litre with Glass Door



Product Code: FPD-09758

Dimensions: 1688mm H x 685mm W x 910mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Interior 12w light - For better visibility of contents without having to open the door
- Front lockable roller castors - To make the fridge easier to move for cleaning and maintenance, and secure it in place
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-R/H only, Replacement Door Seal
- Five wire shelves, Number of Shelf Positions: 30, Size: 536x609x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 40mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR334UK - Pharmacy Fridge 334 Litre with Glass Door



Product Code: FPD-09757

Dimensions: 1860mm H x 595mm W x 660mm D

Internal Dimensions: 1682mm H x 480mm W x 445mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Interior Xw light - For better visibility of contents without having to open the door
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Finish: Painted aluminium
- Colour: White/green
- Door Type: Glass, Hinged-reversible, Door Open Depth: 1180mm, Replacement Door Seal
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Six wire shelves, Number of Shelf Positions: 14, Size: 385x500x6mm
- Automatic defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 220-240v/50hz, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PSR334UK - Pharmacy Fridge 334 Litre with Solid Door



Product Code: FPD-09748

Dimensions: 1860mm H x 595mm W x 670mm D

Internal Dimensions: 1682mm H x 480mm W x 445mm D

- Our large capacity PSR334 solid door refrigerator is ideal for the storage of vaccines and comes complete with our new advanced controller and a fitted lock for added security.
- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Antimicrobial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening.
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Interior Xw light - For better visibility of contents
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Door Type: Solid, Hinged-reversible, Door Open Depth: 1185mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Six wire Shelves, Number of Shelf Positions: 14, Size: 6: 385x500x6mm
- Interior Light (Medical)
- Automatic defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 220-240v/50hz, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PE1607C - Pharmacy Fridge 444 Litre with Solid Door



Product Code: FPD-09749

Dimensions: 1688mm H x 685mm W x 910mm D

- +2°C to +8°C set temperature - Operates at +6°C for optimum storage of temperature controlled drugs and vaccines
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Front lockable roller castors - To make the fridge easier to move for cleaning and maintenance, and secure it in place
- External Finish: Painted steel, Internal Finish: HIPS*
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Type: Solid, Hinged-R/H only, Replacement Door Seal
- Five wire shelves, Number of Shelf Positions: 30, Size: 536x609x6mm
- Automatic defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free, Ammonia Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 40mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM360G - 360 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-04356

Dimensions: 1595mm H x 600mm W x 600mm D

Internal Dimensions: 1418mm H x 513mm W x 441mm D

- Glass Door
- Interior Light
- Gross Volume: 352 Litres, Net Volume 345 Litres
- 5 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48
HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP_™

Shoreline SM360 - 360 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-04364

Dimensions: 1595mm H x 600mm W x 600mm D

Internal Dimensions: 1418mm H x 513mm W x 441mm D

- Capacity 360 Litres Gross, 352 Litres Net
- 5 x Shelves
- 220-230V/50Hz
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- Automatic Defrosting
- Door Lock and keys
- Meets W.H.O. performance criteria (E3/PROC/3/2)
- CFC free
- Reversible Doors
- Adjustable Shelving
- 2-Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM541G - 544 Litres Pharmacy Fridge with Glass Door



Product Code: FPD-07643

Dimensions: 1745mm H x 750mm W x 700mm D

Internal Dimensions: 1452mm H x 600mm W x 560mm D

- High Efficiency Model
- 544 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Remote Alarm Terminals
- Automatic Defrosting
- Door Lock and 2 keys
- Castors (2 front locking)
- Self-Closing Door
- High Efficiency Model
- Thick Wall Insulation
- 1.4 kWhr/24hr
- Adjustable Shelving x 4
- Interior Light
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +38°C
- Made in England
- 2 Years Part & Labour Warranty (UK only)
- Price Includes Delivery to Mainland UK

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Labcold IntelliCold™ RLDG1510 - 430 Litre Pharmacy Fridge with Glass Door



Product Code: FPD-09844

Dimensions: 1860mm H x 600mm W x 700mm D

- A 430 litre glass door pharmacy refrigerator powered by the Labcold IntelliCold™ controller. Two internal fans are fitted to this refrigerator, giving it superior temperature uniformity making it one of the safest ways to store large quantities of cold chain pharmaceuticals. The large safety glass door reduces the need for door openings as the internal light can be operated from the exterior of the fridge so can be turned on independently of the door opening. The easy four button operation, plus the capacity to store temperature data for up to 13 months assist with regulatory compliance and for extra peace of mind it is backed by the Labcold FREE 5 year parts and labour warranty.
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Upright
- Door Type: Glass
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- FREE SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Six Shelves
- Access port



Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold IntelliCold™ RLDF1510 - 430 Litre Pharmacy Fridge with Solid Door



Product Code: FPD-09843

Dimensions: 1860mm H x 600mm W x 700mm D

- A slim line 430 litre pharmacy refrigerator powered by the Labcold IntelliCold™ controller and fitted with two fans for the optimum in temperature stability. Designed exclusively for the safe storage of vaccines and other temperature sensitive pharmaceuticals, this fridge makes compliance with regulatory guidelines simple by featuring an intuitive four button controller and 13 month temperature data logging. A free SD card for temperature data downloads is available when you register for your FREE 5 year parts and labour warranty!
- Can be fitted with Labcold divided drawer(s) which are exclusively designed to keep air flow to a maximum whilst helping to keep stock tidy. Manufactured from hygienic aluminium these drawers can be bought individually and can be retro fitted.
- Configuration: Upright
- Door Type: Solid
- Operating temp (nominal): 5°C
- Temperature range: 2 to 8°C
- Fitted with an IntelliCold™ controller
- 24/7 automatic temperature data logging
- Controller and alarm battery back up door open alarm
- Power failure alarm, High temperature alarm, Low temperature alarm
- Minimum/maximum temperature recording and display
- Off cycle auto defrost
- Fan assisted cooling
- Environmentally friendly refrigerant
- Ammonia free
- Door lock with two keys
- Real time temperature display
- Internal light
- Free SD card when warranty registered
- FREE five year parts and labour warranty
- Compliant with all leading professional guidelines
- Low Energy
- Six Shelves
- Access port



Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SM541 - 544 Litres Pharmacy Fridge with Solid Door



Product Code: FPD-07642

Dimensions: 1745mm H x 750mm W x 700mm D

Internal Dimensions: 1452mm H x 600mm W x 560mm D

- High Efficiency Model
- 544 litres Capacity
- Temperature Range +2°C to +8°C
- Suitable for Vaccine & Pharmaceutical Storage
- Forced-Air Cooling for even temperature distribution
- Digital LED Temperature Display
- Audio / Visual Alarm on Temperature Deviation
- MAX/MIN Temperature Memory
- Remote Alarm Terminals
- Automatic Defrosting
- Door Lock and 2 keys
- Castors (2 front locking)
- Self-Closing Door
- High Efficiency Model
- Thick Wall Insulation
- 1.0 kWhr/24hr
- Adjustable Shelving x 5
- CFC & Ammonia free
- RPSGB, NHS, WHO, RCVS compliant
- 220-230v 50Hz
- Max. Ambient +43°C
- Made in England
- 2 Years Part & Labour Warranty (UK only)
- Price Includes Delivery to Mainland UK

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

LEC PSR600UK - Pharmacy Fridge 600 Litre with Solid Door



Product Code: FPD-10199

Dimensions: 1995mm H x 690mm W x 840mm D

Internal Dimensions: 1386mm H x 560mm W x 653mm D

- Conforms to MHRA guidelines.
- NEW Advanced Controller
- Records temperatures accurate to 0.1°C
- Fitted lock as standard maintains security and ensures the unit can only be accessed by approved personnel
- Fitted Castors for ease of movement
- Wire shelves and basket enables optimum air flow inside the unit
- Temperature range +2°C to +8°C
- Colour White/green
- External finish White
- Internal finish Stainless steel
- Type of plastic SRL800 (HIPS/ABS Mix)
- Door finish Painted stainless steel
- Door type Solid, Hinged
- Door open depth 1145mm
- Adjustable feet (number & position) Fitted castors
- Temperature display on panel
- Temperature alarm
- Five shelves, 22 available positions
- Shelf size 645mm x 523mm
- Shelf material Plastic coated wire
- One plastic coated wire baskets
- Fan
- Automatic defrost
- External digital controller
- Climate class N
- Refrigerant type R134a
- CFC, HFC, HCFC & Ammonia free
- UK 3 pin plug, Other plug types available, Voltage supply 230V, Aux output
- Insulation thickness: 60mm, Insulation type: Cyclopentane
- Optional drawer system available

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR600UK - Pharmacy Fridge 600 Litre with Glass Door



Product Code: FPD-10201

Dimensions: 1995mm H x 690mm W x 840mm D

Internal Dimensions: 1386mm H x 560mm W x 653mm D

- Conforms to MHRA guidelines
- NEW Advanced Controller
- Records temperatures accurate to 0.1°C
- Fitted lock as standard maintains security and ensures the unit can only be accessed by approved personnel
- Fitted castors for ease of movement
- Wire Shelves and Basket enables optimum air flow inside the unit
- Temperature range +2°C to +8°C
- Colour White/green
- External finish White
- Internal finish Stainless steel
- Type of plastic SRL800 (HIPS/ABS Mix)
- Door finish Painted stainless steel
- Door type Glass, Hinged
- Door open depth 1145mm
- Adjustable feet (number & position) Fitted castors
- Temperature display on panel
- Temperature alarm
- Five shelves with 22 available positions
- Shelf size 645mm x 523mm
- One plastic coated wire basket
- Interior light
- Fan
- Automatic defrost
- External digital controller
- Climate class N
- Refrigerant type R134a
- CFC, HFC, HCFC & Ammonia free
- UK 3 pin plug, Other plug types available, Voltage supply 230V, Aux output
- Insulation thickness 60mm, Insulation type Cyclopentane
- Optional drawer system available

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PGR1200UK - Pharmacy Fridge 1200 Litre with Glass Door



Product Code: FPD-10202

Dimensions: 1995mm H x 1330mm W x 840mm D

Internal Dimensions: 1386mm H x 1220mm W x 653mm D

- Conforms to MHRA guidelines
- NEW Advanced Controller records temperatures accurate to 0.1°C
- Fitted Lock as Standard maintains security and ensures the unit can only be accessed by approved personnel
- Fitted Castors for ease of movement
- Wire Shelves and Basket enables optimum air flow inside the unit
- Temperature range +2°C to +8°C
- Colour White/green
- External finish White
- Internal finish Stainless steel
- Type of plastic SRL800 (HIPS/ABS Mix)
- Door finish Painted stainless steel
- Door type Glass, Hinged
- Door open depth 1145mm
- Adjustable feet (number & position) Fitted castors
- Temperature display on panel
- Temperature alarm
- Ten shelves 44 positions available
- Shelf size 645mm x 523mm
- Two plastic coated wire baskets
- Interior light
- Fan
- Automatic defrost
- External digital controller
- Climate class N
- Refrigerant type R134a
- CFC free, HFC, HCFC & Ammonia free
- UK 3 pin plug, Other plug types available, Voltage supply 230V, Aux output
- Insulation thickness 60mm, Insulation type Cyclopentane
- Optional drawer system available

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC PSR1200UK - Pharmacy Fridge 1200 Litre with Solid Door



Product Code: FPD-10200

Dimensions: 1995mm H x 1330mm W x 840mm D

Internal Dimensions: 1386mm H x 1220mm W x 653mm D

- Conforms to MHRA guidelines
- NEW Advanced Controller
- Records temperatures accurate to 0.1°C
- Fitted lock as standard maintains security and ensures the unit can only be accessed by approved personnel
- Fitted castors for ease of movement
- Wire shelves and basket enables optimum air flow inside the unit
- Temperature range +2°C to +8°C
- Colour White/green
- External finish White
- Internal finish Stainless steel
- Type of plastic SRL800 (HIPS/ABS Mix)
- Door finish Painted stainless steel
- Door type Solid, Hinged
- Door open depth 1145mm
- Adjustable feet (number & position) Fitted castors
- Temperature display on panel
- Temperature alarm
- Ten shelves, 44 positions available
- Shelf size 645mm x 523mm
- Two plastic coated wire baskets
- Fan
- Automatic defrost
- External digital controller
- Climate class N
- Refrigerant type R134a
- CFC, HFC, HCFC & Ammonia
- UK 3 pin plug, Other plug types available, Voltage supply 230V, Aux output
- Insulation thickness 60mm, Insulation type Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Advanced Pharmacy Fridges

The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.



[Labcold RPF05043 - 150 Litre
Advanced Pharmacy Fridge with Solid
Door](#)

£1,950.00

Capacity: 150 Litre
Config: Under Bench
820mm H x 600mm W x 620mm D

Product Code: FPD-07551



[Labcold RPF05043 - 150 Litre
Advanced Pharmacy Fridge with Glass
Door](#)

£2,145.00

Capacity: 150 Litre
Config: Under Bench
820mm H x 600mm W x 620mm D

Product Code: FPD-07552



[Labcold RPF21042 - 571 Litre
Advanced Pharmacy Fridge with Solid
Door](#)

£3,555.00

Capacity: 571 Litre
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-07553



[Labcold RPF21042 - 571 Litre
Advanced Pharmacy Fridge with Glass
Door](#)

£4,100.00

Capacity: 571 Litre
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-07554



[Labcold RPF44042 - 1245 Litre
Advanced Pharmacy Fridge with Solid
Door](#)

£5,425.00

Capacity: 1245 Litre
Config: Free Standing
1980mm H x 1440mm W x 820mm D

Product Code: FPD-06566



[Labcold RPF44042 - 1245 Litre
Advanced Pharmacy Fridge with Glass
Door](#)

£6,305.00

Capacity: 1245 Litre
Config: Free Standing
1980mm H x 1440mm W x 820mm D

Product Code: FPD-07555

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF05043 - 150 Litre Advanced Pharmacy Fridge with Solid Door



Product Code: FPD-07551

Dimensions: 820mm H x 600mm W x 620mm D

Internal Dimensions: 680mm H x 540mm W x 350mm D

- Robustly constructed from Stainless Steel this refrigerator from the Labcold Advanced range is designed to provide close temperature control for temperature sensitive products.
- The fan circulated interior has no internal 'step', making it easy to clean and stock and to allow for maximum air flow for even temperature distribution. The two shelves are fully adjustable, so it is possible to maximise the space available.
- Insulated with Envirofoam, which is derived from plant rather than petrochemical oils, the refrigerant is also ammonia, HCFC and CFC free.
- Fitted with a sturdy door lock, this Advanced refrigerator is ideal for busy environments such as busy laboratories or pharmacies.
- Config: Under Counter
- Energy Consumption 3.75 kWh/24h
- Operating Temperature +4°C
- Temperature Set Point Range +2°C to +8°C
- Interior & Exterior Construction Stainless Steel
- Two Shelves
- Shelf dimensions (W xD) 530 x 325mm
- Door Hinging Right hand only
- Door Style: Solid
- Internal Light Glass Door Model only
- Rear Castors
- Door Lock
- Access Port Optional Extra
- Digital Temperature Control, Digital Temperature Display
- High Temperature alarm, Low Temperature alarm, Door Open alarm
- Auto defrost
- Fan Assisted Cooling
- Environmentally friendly insulation
- Eco Friendly Refrigerant
- Ammonia Free
- DH Green Book, RPS/GPhC, NaTHNaC and RCVS Guidelines compliant
- Voltage Supply: 230V/50Hz, UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF21042 - 571 Litre Advanced Pharmacy Fridge with Glass Door



Product Code: FPD-07554

Dimensions: 1980mm H x 720mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 571 Litre
- Auto Defrost
- Shelves: 4
- Operating temp (nominal) +4°C
- Door Type: Glass
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF21042 - 571 Litre Advanced Pharmacy Fridge with Solid Door



Product Code: FPD-07553

Dimensions: 1980mm H x 720mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 571 Litre
- Auto Defrost
- Shelves: 4
- Operating temp (nominal) +4°C
- Door Type: Solid
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF21042 - 571 Litre Advanced Pharmacy Fridge with Glass Door



Product Code: FPD-07554

Dimensions: 1980mm H x 720mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 571 Litre
- Auto Defrost
- Shelves: 4
- Operating temp (nominal) +4°C
- Door Type: Glass
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPFR44042 - 1245 Litre Advanced Pharmacy Fridge with Solid Door



Product Code: FPD-06566

Dimensions: 1980mm H x 1440mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 1245 Litre
- Shelves: 8
- Operating temp (nominal): +4°C
- Door Type: Solid
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPF44042 - 1245 Litre Advanced Pharmacy Fridge with Glass Door



Product Code: FPD-07555

Dimensions: 1980mm H x 1440mm W x 820mm D

- The Labcold Advanced Pharmacy range of refrigerators represents the highest levels of performance, constructional quality, robustness, and appearance available in the market place. Heavy-duty all-metal construction provides strength and durability with ease of cleaning or sanitisation. Microprocessor control, together with ducted airflow and advanced auto defrost system provide class leading temperature stability and uniformity.
- Configuration: Free-standing
- Capacity: 1245 Litre
- Shelves: 8
- Operating temp (nominal): +4°C
- Door Type: Double Glass
- Lockable castors
- Precise microprocessor temperature control
- Digital temperature display and alarm
- Automatic defrost
- Adjustable shelving
- Lockable castors
- Lockable door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Laboratory Refrigeration



Fan Circulated Laboratory Refrigeration
Advanced Laboratory Refrigeration
Sparkfree Laboratory Refrigeration

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Fan Circulated Laboratory Fridges

Fan circulated refrigerators are designed for temperature uniformity throughout the chamber of the refrigerator. Where not integral, all fan circulated fridges can be factory fitted with high/ low temperature alarms and digital temperature displays. Available in a range of sizes, all models feature a door lock for extra security and adjustable shelves so you can make the most of the space available.



[Labcold RLPF05042A Fan Air Circulated Laboratory Fridge 160 Litres](#)

Capacity: 160 Litre
Config: Under Bench
850mm H x 550mm W x 600mm D

Product Code: FPD-07541



[Labcold RLPF09043A Fan Air Circulated Laboratory Fridge 285 Litres](#)

Capacity: 285 Litre
Config: Free Standing
1435mm H x 600mm W x 625mm D

Product Code: FPD-07542



[Labcold RLHF13043 Fan Air Circulated Laboratory Fridge 390 Litres](#)

Capacity: 390 Litre
Config: Free Standing
1880mm H x 580mm W x 615mm D

Product Code: FPD-04922

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPF05042A Fan Air Circulated Laboratory Fridge 160 Litres



Product Code: FPD-07541

Dimensions: 850mm H x 550mm W x 600mm D

- Capacity 160 Litres
- Configuration: Underbench
- Fan assisted air circulation
- Door locks
- CFC, HCFC and Amonia free
- Low energy
- Auto defrosting
- Removable shelves
- Maximum Ambient Temperature up to 32°C
- 3 x Toughened Glass Shelves
- Operating Temperature +4°C
- Temperature Range: +2°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Alarm
- Temperature display

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPF09043A Fan Air Circulated Laboratory Fridge 285 Litres



Product Code: FPD-07542

Dimensions: 1435mm H x 600mm W x 625mm D

- Capacity 285 Litres
- Configuration: Free-standing
- Fan assisted air circulation
- Door locks
- CFC, HCFC and Amonia free
- Low energy
- Auto defrosting
- Removable shelves
- Maximum Ambient Temperature up to 32°C
- 5 Glass Shelves
- Operating Temperature +4°C
- Temperature Range: +2°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Alarm
- Temperature display

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLHF13043 Fan Air Circulated Laboratory Fridge 390 Litres



Product Code: FPD-04922

Dimensions: 1880mm H x 580mm W x 615mm D

- Capacity 390 Litres
- Configuration: Free-standing
- Fan assisted air circulation
- Door locks
- CFC, HCFC and Amonia free
- Temperature Display
- Auto defrosting
- Removable shelves
- Maximum Ambient Temperature up to 32°C
- 6 Shelves
- Operating Temperature +4°C
- Temperature Range: +2°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Alarm
- Temperature display

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Advanced Laboratory Fridges

Advanced fridges are completely constructed from powder coated mild steel with a stainless steel interior so are built to last. Insulated with Envirofoam® and fitted with micro processor controllers these fridges provide accurate temperature control with low GWP. All Advanced refrigerators have optimised air circulation for superior temperature control and feature a temperature display, high and low temperature alarms and door locks fitted as standard. For ease of positioning the larger cabinets have heavy duty wheels with a brake. The Advanced range of refrigerators can be fitted with access ports and remote alarm contacts if required



[Labcold RFR05043 Advanced Laboratory Fridge 150 Litres with Solid Door](#)

£1,905.00

Capacity: 150 Litre
Config: Under Bench
830mm H x 600mm W x 620mm D

Product Code: FPD-10127



[Labcold RFG05043 Advanced Laboratory Fridge 150 Litres with Glass Door](#)

£2,145.00

Capacity: 150 Litre
Config: Under Bench
830mm H x 600mm W x 620mm D

Product Code: FPD-10128



[Labcold RFR21042 Advanced Laboratory Fridge 577 Litres](#)

£3,650.00

Capacity: 577 Litre
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-04944



[Labcold RFG21042 Advanced Laboratory Fridge 577 Litres](#)

£4,210.00

Capacity: 577 Litre
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-04951



[Labcold RFR44042 Advanced Laboratory Fridge 1245 Litres](#)

£5,710.00

Capacity: 1245 Litre
Config: Free Standing
1980mm H x 1440mm W x 820mm D

Product Code: FPD-04945



[Labcold RFG44042 Advanced Laboratory Fridge 1245 Litres](#)

£6,685.00

Capacity: 1245 Litre
Config: Free Standing
1980mm H x 1440mm W x 820mm D

Product Code: FPD-04952

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR05043 Advanced Laboratory Fridge 150 Litres with Solid Door



Product Code: FPD-10127

Dimensions: 830mm H x 600mm W x 620mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Capacity: 150 litres (gross)
- Door Type: Solid
- Construction: AISI 304 grade stainless steel interior and exterior.
- Shelves Supplied: 2, Plastic coated wire mesh 530 x 325mm.
- Microprocessor temperature control with digital display to 1°C
- Hi/Lo temperature alarm - audible
- Fan circulated interior
- Door Lock
- CFC and HFC free
- Temperature set-point range: +2 to +8°C
- Automatic defrost with condensate removal
- Maximum ambient temperature: +43°C

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR05043 Advanced Laboratory Fridge 150 Litres with Solid Door



Product Code: FPD-10127

Dimensions: 830mm H x 600mm W x 620mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Capacity: 150 litres (gross)
- Door Type: Solid
- Construction: AISI 304 grade stainless steel interior and exterior.
- Shelves Supplied: 2, Plastic coated wire mesh 530 x 325mm.
- Microprocessor temperature control with digital display to 1°C
- Hi/Lo temperature alarm - audible
- Fan circulated interior
- Door Lock
- CFC and HFC free
- Temperature set-point range: +2 to +8°C
- Automatic defrost with condensate removal
- Maximum ambient temperature: +43°C

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR21042 Advanced Laboratory Fridge 577 Litres



Product Code: FPD-04944

Dimensions: 1980mm H x 720mm W x 820mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Configuration: Free-standing, solid door
- Capacity: 577 Litres
- Operating temperature: +4°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 4 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RAG21042 Advanced Laboratory Fridge 577 Litres



Product Code: FPD-04951

Dimensions: 1980mm H x 720mm W x 820mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Configuration: Free-standing, solid door
- Capacity: 577 Litres
- Operating temperature: +4°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 4 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RAFG44042 Advanced Laboratory Fridge 1245 Litres



Product Code: FPD-04952

Dimensions: 1980mm H x 1440mm W x 820mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Configuration: Free-standing, double glass door
- Capacity: 1245 Litres
- Operating temperature: +4°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 8 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RAFR44042 Advanced Laboratory Fridge 1245 Litres



Product Code: FPD-04945

Dimensions: 1980mm H x 1440mm W x 820mm D

- For hospital and pharmacy use where close temperature control and uniformity is required together with a high level of construction, robustness and features, and where there is no risk that volatile, flammable materials or components might be stored.
- Configuration: Free-standing, double solid door
- Capacity: 1245 Litres
- Operating temperature: +4°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 8 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Sparkfree Laboratory Fridges

The Sparkfree Fridges is engineered to store flammable materials safely with no sources of ignition inside the chambers. Everything that could conceivably cause a spark, such as the controller, is located outside of the refrigerator chamber making these fridges and freezers an important part of laboratory safety. Designed specifically for laboratories, these refrigerators are solidly built and come in a range of sizes, from a compact bench top to a large 400 litre refrigerator. All Sparkfree fridges are factory fitted with a door lock for optimum security.



[Labcold RLPR02042 Sparkfree Laboratory Fridge 60 Litres](#)

Capacity: 60 Litre
Config: Upright Refrigerator
572mm H x 485mm W x 520mm D
Product Code: FPD-04912



[LEC LR207C Laboratory Refrigerator 82 Litre](#)

Capacity: 82 Litre
Config: Bench Top
660mm H x 502mm W x 540mm D
Product Code: FPD-09794



[LEC LR307C Laboratory Refrigerator 107 Litre](#)

Capacity: 107 Litre
Config: Under Bench
830mm H x 495mm W x 620mm D
Product Code: FPD-09795



[Liebherr LKUexv 1610 Laboratory Refrigerator - 141 Litre](#)

Capacity: 141 Litre
Config: Bench Top
820mm H x 600mm W x 615mm D
Product Code: FPD-10260



[LEC LSR151 Laboratory Refrigerator 151 Litre](#)

Capacity: 151 Litre
Config: Under Bench
845mm H x 595mm W x 595mm D
Product Code: FPD-09796



[Labcold RLPR05042A Sparkfree Laboratory Fridge 160 Litres](#)

Capacity: 160 Litre
Config: Upright Refrigerator
850mm H x 550mm W x 600mm D
Product Code: FPD-07537



[LEC LSR288 Laboratory Refrigerator 288L Litre](#)

Capacity: 284 Litre
Config: Free Standing
1570mm H x 595mm W x 634mm D
Product Code: FPD-09797



[Labcold RLPR09043A Sparkfree Laboratory Fridge 285 Litres](#)

Capacity: 285 Litre
Config: Upright Refrigerator
1435mm H x 600mm W x 625mm D
Product Code: FPD-07538

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Sparkfree Laboratory Fridges

The Sparkfree Fridges is engineered to store flammable materials safely with no sources of ignition inside the chambers. Everything that could conceivably cause a spark, such as the controller, is located outside of the refrigerator chamber making these fridges and freezers an important part of laboratory safety. Designed specifically for laboratories, these refrigerators are solidly built and come in a range of sizes, from a compact bench top to a large 400 litre refrigerator. All Sparkfree fridges are factory fitted with a door lock for optimum security.



[Labcold RLHD13043 Sparkfree Laboratory Fridge 360 Litres](#)

Capacity: 360 Litre
Config: Upright Refrigerator
1750mm H x 585mm W x 655mm D
Product Code: *FPD-04917*



[Labcold RLSP1505 - 400 Litre Sparkfree Laboratory Fridge with Solid Door](#)

Capacity: 400 Litre
Config: Upright Refrigerator
1860mm H x 600mm W x 695mm D
Product Code: *FPD-07707*



[LEC LR1607C Laboratory Refrigerator 474 Litre](#)

Capacity: 474 Litre
Config: Free Standing
1688mm H x 695mm W x 910mm D
Product Code: *FPD-09798*

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPR02042 Sparkfree Laboratory Fridge 60 Litres



Product Code: FPD-04912

Dimensions: 572mm H x 485mm W x 520mm D

- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity 60 Litres
- Configuration: Counter top
- Sparkfree safety interiors
- CFC, HCFC and Amonia free
- Automatic defrosting
- Removable Shelves
- Maximum ambient temperatures up to 32°C
- Door Lock
- 2 Shelves
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temp Display
- Alarm
- Wall Mounting Kit

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LR207C Laboratory Refrigerator 82 Litre



Product Code: FPD-09794

Dimensions: 660mm H x 502mm W x 540mm D

Internal Dimensions: 540mm H x 415mm W x 340mm D

- +2°C to +10°C set temperature - Operates at +6°C for best possible storage of samples and chemicals
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- Small countertop size - For limited space needing biomedical and laboratory refrigeration
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid - Hinged, reversible
- Door Open Depth: 976mm
- Two wire Shelves: 2
- Shelf Size: 438x274x6mm
- Automatic Defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <43dB
- Insulation Thickness: 30mm
- Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 0.77

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LR307C Laboratory Refrigerator 107 Litre



Product Code: FPD-09795

Dimensions: 830mm H x 495mm W x 620mm D

Internal Dimensions: 540mm H x 403mm W x 400mm D

- +2°C to +10°C set temperature - Operates at +6°C for best possible storage of samples and chemicals
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark-free - There is nothing electrical inside the fridge so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Type: Solid - Hinged, reversible
- Door Open Depth: 1050mm
- Two Adjustable Front Feet
- Two Wire Shelves in Four Available Positions
- Shelf Size: 427x374x6mm
- Automatic Defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <43dB
- Insulation Thickness 35mm
- Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 0.86

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Liebherr LGUex 1500 Laboratory Freezer - 139 Litre



Product Code: FPD-10261

Dimensions: 820mm H x 600mm W x 615mm D

- Temperature range -9°C to -26°C
- Suitable for the storage of explosive or highly inflammable substances
- Digital Temperature Control & Display
- Temperature deviation alarm with volt free contacts
- Door open alarm
- Manual Defrost
- Door Lock and 2 keys
- CFC & Ammonia free
- 3 x drawers + 1 Basket
- 220-230v 50Hz
- Max. Ambient +38°C

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Liebherr LKUexv 1610 Laboratory Refrigerator - 141 Litre



Product Code: FPD-10260

Dimensions: 820mm H x 600mm W x 615mm D

- Temperature range +3°C to +8°C
- Suitable for the storage of explosive or highly inflammable substances
- Digital Temperature Control & Display
- Temperature deviation alarm with volt free contacts
- Door open alarm
- Auto Defrost
- Door Lock and 2 keys
- CFC & Ammonia free
- 3 x shelves
- 220-230v 50Hz
- Max. Ambient +38°C
- 2 Year Warranty

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

LEC LSR151 Laboratory Refrigerator 151 Litre



Product Code: FPD-09796

Dimensions: 845mm H x 595mm W x 595mm D

Internal Dimensions: 705mm H x 500mm W x 450mm D

- Set temperature - Operates at +X°C for best possible storage of samples and chemicals
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: +2°C to +10°C
- Colour: white/green
- Door Type: Solid - Hinged, reversible
- Open Door Sensor: Yes
- Lock
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Temperature Display on Panel
- Temperature Alarm
- Three Wire Shelves Six Positions
- Shelf Size: 320x520x6mm
- Two Adjustable Door Wire Bins
- Automatic Defrost
- External Digital Controller
- Climate Class: SN
- Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Insulation Thickness: 45mm
- Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPR05042A Sparkfree Laboratory Fridge 160 Litres



Product Code: FPD-07537

Dimensions: 850mm H x 550mm W x 600mm D

- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity 160 Litres
- Configuration: Underbench lockable
- Sparkfree safety interiors
- CFC, HCFC and Amonia free
- Automatic defrosting
- Removable Shelves
- Maximum ambient temperatures up to 32°C
- Door Lock
- Low Energy
- 4 Toughened Glass Shelves
- 1 x Drawer
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temp Display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSR288 Laboratory Refrigerator 288L Litre



Product Code: FPD-09797

Dimensions: 1570mm H x 595mm W x 634mm D

Internal Dimensions: 1382mm H x 480mm W x 440mm D

- Set temperature - Operates at +X°C for best possible storage of samples and chemicals
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- External Finish: white
- Internal Finish: white
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: +2°C to +10°C
- Colour: White/green
- Door Type: Solid - Hinged, reversible
- Door Open Depth: 1185mm
- Open Door Sensor
- Lock
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Temperature Display on Panel
- Temperature Alarm
- Five Wire Shelves, 12 Positions
- Shelf Size: 540x497x6mm
- One Solid Moulded Drawer
- Four Solid Adjustable Door Bins
- Automatic Defrost
- External Digital Controller
- Climate Class: SN
- Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Insulation Thickness: 60mm
- Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLPR09043A Sparkfree Laboratory Fridge 285 Litres



Product Code: FPD-07538

Dimensions: 1435mm H x 600mm W x 625mm D

- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity 285 Litres
- Configuration: Free-standing lockable
- Sparkfree safety interiors
- CFC, HCFC and Amonia free
- Automatic defrosting
- Removable Shelves
- Maximum ambient temperatures up to 32°C
- Low Energy
- Door Lock
- 5 Glass Shelves
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temp Display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLHD13043 Sparkfree Laboratory Fridge 360 Litres



Product Code: FPD-04917

Dimensions: 1750mm H x 585mm W x 655mm D

- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity 360 Litres
- Configuration: Free-standing
- Sparkfree safety interiors
- CFC, HCFC and Amonia free
- Automatic defrosting
- Removable Shelves
- Maximum ambient temperatures up to 32°C
- 6 Shelves
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temp Display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLSP1505 - 400 Litre Sparkfree Laboratory Fridge with Solid Door



Product Code: FPD-07707

Dimensions: 1860mm H x 600mm W x 695mm D

- A purpose-built refrigerator for laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future. Fan circulation of chamber air provides much improved temperature uniformity.
- Capacity: 400 litres (gross)
- Door Type: Solid
- Construction: White coated steel interior and exterior
- Shelves Supplied: 6, Plastic coated wire mesh
- Sparkfree Interior
- Fan assisted air circulation powered by external motor
- Digital display of temperature
- Audible alarm with mute facility
- Operating temperature +4°C
- Temperature range: 0°C to +10°C
- Door Lock
- Easy-grip door handle
- Access Port
- CFC, HCFC and Ammonia free
- Accurate temperature
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LR1607C Laboratory Refrigerator 474 Litre



Product Code: FPD-09798

Dimensions: 1688mm H x 695mm W x 910mm D

Internal Dimensions: 1170mm H x 550mm W x 670mm D

- +2°C to +10°C set temperature - Operates at +6°C for best possible storage of samples and chemicals
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Front lockable roller castors - To make this large fridge easier to move for cleaning and maintenance
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Type: Solid - Hinged, R/H only
- Door Open Depth: 1596mm
- Lock
- Adjustable Feet: Castors
- Five Wire Shelves 30 Positions: 30
- Shelf Size: 536x609x6mm
- Automatic Defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Insulation Thickness: 65mm
- Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Laboratory Freezers



Sparkfree Laboratory Freezers
Advanced Laboratory Freezers
Sparkfree Laboratory Fridge/Freezers

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Sparkfree Laboratory Freezers

The Labcold Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Everything that could conceivably cause a spark, such as the controller, is located outside of the freezer chamber making these freezers an important part of laboratory safety. Designed specifically for laboratories, these freezers are solidly built and come in both upright and chest configuration. All Labcold Sparkfree freezers and chest freezers are factory fitted with a door lock for optimum security.



[Labcold RLVF02203 Sparkfree Laboratory Freezer 36 Litres](#)

Capacity: 36 Litre
Config: Upright Freezer
515mm H x 440mm W x 500mm D
Product Code: FPD-09917



[Labcold RLVF02203/LK Sparkfree Laboratory Freezer 36 Litres Lockable](#)

Capacity: 36 Litre
Config: Upright Freezer
515mm H x 440mm W x 500mm D
Product Code: FPD-09918



[LEC ISU27C Laboratory Freezer 50 Litres](#)

Capacity: 50 Litre
Config: Bench Top
562mm H x 500mm W x 510mm D
Product Code: FPD-09804



[LEC ISU37C Laboratory Freezer 102 Litres](#)

Capacity: 102 Litre
Config: Under Bench
830mm H x 495mm W x 620mm D
Product Code: FPD-09805



[Labcold RLVF04202A Sparkfree Laboratory Freezer 120 Litres](#)

Capacity: 120 Litre
Config: Upright Freezer
850mm H x 545mm W x 565mm D
Product Code: FPD-07539



[Labcold RLVF04202A/LK Sparkfree Laboratory Freezer 120 Litres](#)

Capacity: 120 Litre
Config: Upright Freezer
850mm H x 545mm W x 565mm D
Product Code: FPD-07540



[Liebherr LGUex 1500 Laboratory Freezer - 139 Litre](#)

Capacity: 139 Litre
Config: Bench Top
820mm H x 600mm W x 615mm D
Product Code: FPD-10261



[LEC LSF151UK Laboratory Freezer 151 Litres](#)

Capacity: 151 Litre
Config: Under Bench
845mm H x 595mm W x 595mm D
Product Code: FPD-09806

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Sparkfree Laboratory Freezers

The Labcold Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Everything that could conceivably cause a spark, such as the controller, is located outside of the freezer chamber making these freezers an important part of laboratory safety. Designed specifically for laboratories, these freezers are solidly built and come in both upright and chest configuration. All Labcold Sparkfree freezers and chest freezers are factory fitted with a door lock for optimum security.



[Labcold RLVF07203 - 180 Litre Sparkfree Laboratory Freezer with Solid Door](#)

Capacity: 180 Litre
Config: Upright Freezer
1440mm H x 544mm W x 570mm D
Product Code: FPD-07708



[Labcold RLVF07203LK - 180 Litre Sparkfree Laboratory Freezer with Solid Door](#)

Capacity: 180 Litre
Config: Upright Freezer
1440mm H x 544mm W x 570mm D
Product Code: FPD-07709



[Labcold RLCF0720 Sparkfree Laboratory Freezer 215 Litres](#)

Capacity: 215 Litre
Config: Chest Freezer
860mm H x 720mm W x 730mm D
Product Code: FPD-04934



[LEC LSF232UK Laboratory Freezer 232 Litres](#)

Capacity: 232 Litre
Config: Free Standing
1570mm H x 595mm W x 634mm D
Product Code: FPD-09807



[Labcold RLVF14201 Sparkfree Laboratory Freezer 375 Litres](#)

Capacity: 375 Litre
Config: Upright Freezer
1850mm H x 595mm W x 640mm D
Product Code: FPD-04933



[Labcold RLCF1520 Sparkfree Laboratory Freezer 447 Litres](#)

Capacity: 447 Litre
Config: Chest Freezer
860mm H x 1300mm W x 730mm D
Product Code: FPD-04935



[Labcold RLCF2120 Sparkfree Laboratory Freezer 607 Litres](#)

Capacity: 607 Litre
Config: Chest Freezer
860mm H x 1700mm W x 730mm D
Product Code: FPD-04936

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF02203 Sparkfree Laboratory Freezer 36 Litres

Product Code: FPD-09917

Dimensions: 515mm H x 440mm W x 500mm D



- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future
- Capacity: 36 Litres (gross)
- Construction: White enamel exterior, fully moulded plastic interior
- Shelves supplied: 1, wire mesh
- Sparkfree interior
- Counter- top use
- CFC and HCFC free
- Temperature set-point range: -10 to -26oC
- Operator controlled defrost

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ISU27C Laboratory Freezer 50 Litres



Product Code: FPD-09804

Dimensions: 562mm H x 500mm W x 510mm D

Internal Dimensions: 365mm H x 345mm W x 330mm D

- -18°C to -20°C set temperature - Perfect for short term sample storage when space is at a premium
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark-free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- Small countertop size - For limited space needing biomedical and laboratory refrigeration
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid
- Door Style: Hinged, reversible
- Door Open Depth: 976mm
- Lock: Yes
- One fixed wire shelf
- Manual defrost
- Climate Class: N
- Refrigerant Type: R600a
- CFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <45dB
- Insulation Thickness: 30mm
- Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ISU37C Laboratory Freezer 102 Litres



Product Code: FPD-09805

Dimensions: 830mm H x 495mm W x 620mm D

Internal Dimensions: 683mm H x 414mm W x 415mm D

- -18°C to -20°C set temperature - Perfect for short term sample storage when space is at a premium
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid, Hinged-reversible
- Door Open Depth: 1050mm
- Lock: Yes
- Adjustable Feet: 1, front
- Two adjustable wire shelves
- Manual defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free & HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <43dB
- Insulation Thickness: 35mm
- Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.3

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF04202A/LK Sparkfree Laboratory Freezer 120 Litres



Product Code: FPD-07540

Dimensions: 850mm H x 545mm W x 565mm D

- Capacity 120 Litres
- Configuration: Under-bench lockable
- Sparkfree safety interior
- CFC, HCFC and Amonia free
- Operating temperature -20°C
- Operator controlled defrost
- Maximum ambient temperatures up to 32°C
- Low energy
- Door lock
- 4 Shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temperature display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF04202A Sparkfree Laboratory Freezer 120 Litres



Product Code: FPD-07539

Dimensions: 850mm H x 545mm W x 565mm D

- Capacity 120 Litres
- Configuration: Under-bench
- Sparkfree safety interior
- CFC, HCFC and Amonia free
- Operating temperature -20°C
- Operator controlled defrost
- Maximum ambient temperatures up to 32°C
- Low energy
- 4 Plastic removeable drawers separated by fixed refrigerated shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Temperature display
- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Liebherr LGUex 1500 Laboratory Freezer - 139 Litre



Product Code: FPD-10261

Dimensions: 820mm H x 600mm W x 615mm D

- Temperature range -9°C to -26°C
- Suitable for the storage of explosive or highly inflammable substances
- Digital Temperature Control & Display
- Temperature deviation alarm with volt free contacts
- Door open alarm
- Manual Defrost
- Door Lock and 2 keys
- CFC & Ammonia free
- 3 x drawers + 1 Basket
- 220-230v 50Hz
- Max. Ambient +38°C

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

LEC LSF151UK Laboratory Freezer 151 Litres



Product Code: FPD-09806

Dimensions: 845mm H x 595mm W x 595mm D

- -18°C to -25°C set temperature - Perfect for short term sample storage when space it at a premium
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Roller castors - To make the freezer easier to move for cleaning and maintenance
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Open Depth: 1165mm
- Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three fixed wire basket Shelves: 3
- Manual defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Aux Output
- Insulation Thickness: 45mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF07203LK - 180 Litre Sparkfree Laboratory Freezer with Solid Door



Product Code: FPD-07709

Dimensions: 1440mm H x 544mm W x 570mm D

- Capacity: 180 Litres (gross)
- Construction: White enamel exterior, plastic interior
- Door Type: Solid
- Internal Arrangement: 4 drawers and 2 compartments with pull down flaps separated by fixed refrigerated shelves
- Sparkfree Interior
- CFC and HCFC free
- Temperatures set-point range: -18oC to -26oC
- Analogue thermostat control
- Visual alarm for high temperature
- Low energy operation
- Operator controlled defrost
- Locking Door
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF07203 - 180 Litre Sparkfree Laboratory Freezer with Solid Door



Product Code: FPD-07708

Dimensions: 1440mm H x 544mm W x 570mm D

- Capacity: 180 Litres (gross)
- Construction: White enamel exterior, plastic interior
- Door Type: Solid
- Internal Arrangement: 4 drawers and 2 compartments with pull down flaps separated by fixed refrigerated shelves
- Sparkfree Interior
- CFC and HCFC free
- Temperatures set-point range: -18oC to -26oC
- Analogue thermostat control
- Visual alarm for high temperature
- Low energy operation
- Operator controlled defrost
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLCF0720 Sparkfree Laboratory Freezer 215 Litres



Product Code: FPD-04934

Dimensions: 860mm H x 720mm W x 730mm D

- Configuration: Chest
- Capacity 215 Litres
- Operating temperature -20°C
- Sparkfree interior
- Door lock
- Alarm
- Low energy
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSF232UK Laboratory Freezer 232 Litres



Product Code: FPD-09807

Dimensions: 1570mm H x 595mm W x 634mm D

- -18°C to -25°C Set temperature - Perfect for short term storage for a large amount of samples
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark-free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Roller castors - To make the freezer easier to move for cleaning and maintenance
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Colour: White/green
- Door Open Depth: 1185mm
- Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Four wire drawers, 2 fixed flaps, Four Moulded Drawer
- Manual defrost
- Climate Class: SN, Refrigerant Type: R600a
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Aux Output: Yes
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF14201 Sparkfree Laboratory Freezer 375 Litres



Product Code: FPD-04933

Dimensions: 1850mm H x 595mm W x 640mm D

- Capacity 375 Litres
- Configuration: Free-standing
- Sparkfree safety interior
- CFC, HCFC and Amonia free
- Operating temperature -20°C
- Operator controlled defrost
- Maximum ambient temperatures up to 32°C
- 6 Shelves
- Temperature Display
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Optional Extras

- Alarm

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLCF1520 Sparkfree Laboratory Freezer 447 Litres



Product Code: FPD-04935

Dimensions: 860mm H x 1300mm W x 730mm D

- Configuration: Chest
- Capacity 447 Litres
- Operating temperature -20°C
- Sparkfree interior
- Door lock
- Alarm
- Low energy
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLCF2120 Sparkfree Laboratory Freezer 607 Litres



Product Code: FPD-04936

Dimensions: 860mm H x 1700mm W x 730mm D

- Configuration: Chest
- Capacity 607 Litres
- Operating temperature -20°C
- Sparkfree interior
- Door lock
- Alarm
- Low energy
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Advanced Laboratory Freezers

Advanced freezers are completely constructed from powder coated mild steel with a stainless steel interior so are built to last. Insulated with Envirofoam and fitted with micro processor controllers these freezers provide accurate temperature control with low GWP. All Advanced freezers have optimised air circulation for superior temperature control and feature a temperature display, high and low temperature alarms and door locks fitted as standard. For ease of positioning the larger cabinets have heavy duty wheels with a brake. Designed for any situation where accurate temperature control is vital and there is a need to store a large quantity of product. Advanced freezers can be fitted with access ports and remote alarm contacts.




[Labcold RAFR05202 Advanced Laboratory Freezer 150 Litres](#)

Capacity: 150 Litre
Config: Upright Freezer
830mm H x 600mm W x 620mm D
Product Code: FPD-04946



[Labcold RAFR21262 Advanced Laboratory Freezer 571 Litres](#)

Capacity: 571 Litre
Config: Upright Freezer
1980mm H x 720mm W x 820mm D
Product Code: FPD-04948



[Labcold RAFR44262 Advanced Laboratory Freezer 1245 Litres](#)

Capacity: 1245 Litre
Config: Upright Freezer
1980mm H x 1440mm W x 820mm D
Product Code: FPD-04949

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR05202 Advanced Laboratory Freezer 150 Litres



Product Code: FPD-04946

Dimensions: 830mm H x 600mm W x 620mm D

- Configuration: Underbench, solid door
- Capacity: 150 Litres
- Operating temperature: -20°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 3 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RAFR44262 Advanced Laboratory Freezer 1245 Litres



Product Code: FPD-04949

Dimensions: 1980mm H x 1440mm W x 820mm D

- Configuration: Free-standing, double solid door
- Capacity: 1245 Litres
- Operating temperature: -20°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 8 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RFR21262 Advanced Laboratory Freezer 571 Litres



Product Code: FPD-04948

Dimensions: 1980mm H x 720mm W x 820mm D

- Configuration: Free-standing, solid door
- Capacity: 571 Litres
- Operating temperature: -20°C
- Precise microprocessor temperature control
- Digital temperature display
- Audible and visual temperature alarm
- Optimised air circulation and temperature uniformity
- Automatic defrosting with minimal temperature deviation
- Lockable castors
- Heavy duty all-metal construction
- Door lock
- 4 Adjustable heavy duty shelves
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Combined Laboratory Fridge/Freezers

This is the latest generation of combined sparkfree laboratory refrigerator and freezer. Providing all the routine storage required on a floor area no bigger than a standard fridge. For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future. Both refrigerated and frozen samples can be stored and accessed within the same small floor area.



[LEC ISR27C Laboratory Fridge Freezer](#)

Capacity: 83 Litre
Config: Bench Top
660mm H x 502mm W x 540mm D
Product Code: FPD-09808



[LEC LSC119UK Laboratory Fridge Freezer](#)

Capacity: 119 Litre
Config: Under Bench
860mm H x 540mm W x 600mm D
Product Code: FPD-09809



[LEC LSC263UK Laboratory Fridge Freezer](#)

Capacity: 263 Litre
Config: Upright Refrigerator
1680mm H x 540mm W x 600mm D
Product Code: FPD-09810



[LEC LSC324UK Laboratory Fridge Freezer](#)

Capacity: 324 Litre
Config: Upright Refrigerator
1870mm H x 600mm W x 600mm D
Product Code: FPD-09811



[Labcold RLFF13246 Laboratory Fridge/Freezer 230/160 Litres](#)

Capacity: 390 Litre
1965mm H x 595mm W x 630mm D
Product Code: FPD-04937



[Labcold RLFF13246/LK Laboratory Fridge/Freezer 230/160 Litres](#)

Capacity: 390 Litre
1965mm H x 595mm W x 630mm D
Product Code: FPD-04938

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLVF02203/LK Sparkfree Laboratory Freezer 36 Litres Lockable

Product Code: FPD-09918

Dimensions: 515mm H x 440mm W x 500mm D



- For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future
- Capacity: 36 Litres (gross)
- Construction: White enamel exterior, fully moulded plastic interior
- Shelves supplied: 1, wire mesh
- Sparkfree interior
- Counter- top use
- CFC and HCFC free
- Temperature set-point range: -10 to -26oC
- Operator controlled defrost

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ISR27C Laboratory Fridge Freezer



Product Code: FPD-09808

Dimensions: 660mm H x 502mm W x 540mm D

- +2°C to +10°C fridge set temperature, -6°C to -9°C freezer set temperature - Designed for when both refrigeration and freezing is needed, perfect for smaller spaces
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge and freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge or freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- Small countertop size - For limited space needing biomedical and laboratory refrigeration
- Fridge Net Capacity: 75, Freezer Net Capacity: 5.5
- External Finish: White
- Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS Mix)
- Colour: White/green
- Door Type: Solid
- Door Open Depth: 976mm
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Two wire Shelves, Shelf Size: 438x274x6mm & 438x137x6mm
- Fridge: Automatic defrost, Freezer: Manual defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, CFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Spark-free
- Noise Level: <43dB
- Insulation Thickness: 30mm, Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.74

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSC119UK Laboratory Fridge Freezer



Product Code: FPD-09809

Dimensions: 860mm H x 540mm W x 600mm D

- Set temperature - Designed for when both refrigeration and freezing is needed, perfect for smaller spaces
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge and freezer
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge-freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Roller castors - To make the fridge-freezer easier to move for cleaning and maintenance
- Fridge Net Capacity: 105L, Freezer Net Capacity: 12L
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: Fridge +2°C to +10°C, freezer -18°C to -25°C
- Colour: White/green
- Door Type: Solid, Door Open Depth: 1090mm, Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Two wire shelves & One solid moulded Drawer, Two full adjustable door bins + two half
- Fridge: Automatic defrost, Freezer: Manual defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Aux Output: Yes
- Insulation Thickness: 45mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSC263UK Laboratory Fridge Freezer



Product Code: FPD-09810

Dimensions: 1680mm H x 540mm W x 600mm D

- Our spark free laboratory fridge-freezer offers great flexibility and is ideal for the safe storage of a wide variety of volatile liquids and chemicals.
- Set temperature - Designed for when both refrigeration and freezing is needed, perfect for smaller spaces
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge and freezer
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the fridge-freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Reversible door - For maximum versatility in sometimes limited working spaces
- Roller castors - To make the fridge-freezer easier to move for cleaning and maintenance
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: Fridge +2°C to +10°C, freezer -18°C to -25°C
- Door Type: Solid, Door Open Depth: 1110mm, Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three wire Shelves: 3 (fridge), Moulded Drawer: Fridge: 2 solid Adjustable Door Bins: 2 full + 2 half
- Fridge: Automatic defrost, Freezer: Manual defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane
- Fridge Net Capacity: 157 Litre, Freezer Net Capacity: 88 Litre

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC LSC324UK Laboratory Fridge Freezer



Product Code: FPD-09811

Dimensions: 1870mm H x 600mm W x 600mm D

- Our new LSC324UK is our largest fridge-freezer unit yet, offering you even greater storage capacity
- Set temperature - Designed for when both refrigeration and freezing is needed, perfect for smaller spaces
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge and freezer
- Antimicrobial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark-free - There is nothing electrical inside the fridge-freezer, so it is safe to store reactive substances in without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Reversible door - for maximum versatility in sometimes limited working spaces
- Roller castors - To make the fridge-freezer easier to move for cleaning and maintenance
- External Finish: White, Internal Finish: White
- Type of plastic: SRL800 (HIPS / ABS mix)
- Temperature Range: Fridge +2°C to +10°C, freezer -18°C to -25°C
- Door Type: Solid, Door Open Depth: 1190mm, Open Door Sensor
- Adjustable Feet: 2 adjustable fronts, 2 rear rollers
- Three wire Shelves, Two solid moulded drawers, Adjustable Door Bins: 2 full + 4 half
- Fridge: Automatic defrost, Freezer: Manual defrost
- Climate Class: SN
- CFC Free, HFC Free, HCFC
- UK 3 Pin Plug, Voltage Supply: 230V, Aux Output
- Insulation Thickness: 60mm, Insulation Type: Cyclopentane

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLFF13246/LK Laboratory Fridge/Freezer 230/160 Litres



Product Code: FPD-04938

Dimensions: 1965mm H x 595mm W x 630mm D

- Lab Fridge Freezer
- Configuration: Freestanding
- Capacity 170/160 Litres
- +4/-20°C
- Door lock
- Sparkfree interior
- Auto defrosting on freezer
- Low energy
- 4 shelves in refrigeration
- 5 shelves in freezer
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLFF13246 Laboratory Fridge/Freezer 230/160 Litres



Product Code: FPD-04937

Dimensions: 1965mm H x 595mm W x 630mm D

- Lab Fridge Freezer
- Configuration: Freestanding
- Capacity 170/160 Litres
- +4/-20°C
- Sparkfree interior
- Auto defrosting on freezer
- Low energy
- 4 shelves in refrigeration
- 5 shelves in freezer
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696

Fax 0115 727 0123

Email info@pharmacy-fridge.co.uk

Low Temperature Refrigeration



Low Temperature Chest Freezers
Low Temperature Freezers
Ultra Low Temperature Freezers
Laboratory Super Freezers
Compact Low Temperature Laboratory Freezers
DuoCold Ultra Low Temperature Freezers

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Low Temperature Chest Freezers

A variety of specialist Low Temperature Chest Freezers Minus 10°C to Minus 45°C for the storage of chemicals, samples and other important matter.



[LEC CLT140 Low Temperature 140 Litre Chest Freezer](#)

Capacity: 140 Litre
Config: Chest
885mm H x 720mm W x 600mm D

Product Code: FPD-10214



[LEC CLT296 Low Temperature 296 Litre Chest Freezer](#)

Capacity: 296 Litre
Config: Chest
885mm H x 1260mm W x 600mm D

Product Code: FPD-10215



[LEC CLT383 Low Temperature 383 Litre Chest Freezer](#)

Capacity: 383 Litre
Config: Chest
885mm H x 1560mm W x 600mm D

Product Code: FPD-10216

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC CLT140 Low Temperature 140 Litre Chest Freezer



Product Code: FPD-10214

Dimensions: 885mm H x 720mm W x 600mm D

Internal Dimensions: 650mm H x 560mm W x 440mm D

- -10°C to -45°C set temperature - Ideal for medium term sample storage
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Colour: White
- Door Type: solid
- Door Style: hinged lid
- Lock
- Temperature Display on Panel
- Temperature Alarm
- Manual defrost: manual
- Controller type: digital, external
- Climate Class: N
- Refrigerant Type: R414
- CFC Free
- UK 3 Pin Plug
- Voltage Supply: 230/240V
- Spark-free
- Aux Output

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC CLT296 Low Temperature 296 Litre Chest Freezer



Product Code: FPD-10215

Dimensions: 885mm H x 1260mm W x 600mm D

Internal Dimensions: 650mm H x 1100mm W x 440mm D

- -10° to -45° set temperature - Ideal for medium term sample storage
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Colour: White
- Door Type: solid
- Door Style: hinged lid
- Lock
- Temperature Display on Panel
- Temperature Alarm
- Auto or Manual defrost: manual
- Controller type: digital, external
- Climate Class: N
- Refrigerant Type: R414
- CFC Free
- UK 3 Pin Plug
- Voltage Supply: 230/240
- Spark-free
- Aux Output

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC CLT383 Low Temperature 383 Litre Chest Freezer



Product Code: FPD-10216

Dimensions: 885mm H x 1560mm W x 600mm D

Internal Dimensions: 650mm H x 1400mm W x 440mm D

- -10°C to -45°C set temperature - Ideal for medium term sample storage
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Colour: White
- Door Type: solid
- Door Style: hinged lid
- Lock
- Temperature Display on Panel
- Temperature Alarm
- Manual defrost
- Controller type: digital, external
- Climate Class: N
- Refrigerant Type: R414
- CFC Free
- UK 3 Pin Plug
- Voltage Supply: 230/240V
- Spark-free
- Aux Output

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Low Temperature Freezers

A variety of specialist Low Temperature Freezers Minus 20°C to Minus 40°C for the storage of chemicals, samples and other important matter.



[LEC CLT111 Low Temperature 113 Litre Chest Freezer](#)

Capacity: 113 Litre
Config: Chest
1080mm H x 730mm W x 850mm D
Product Code: FPD-10217



[LEC ULT341 Low Temperature 322 Litre Upright Freezer](#)

Capacity: 322 Litre
Config: Free Standing
1985mm H x 765mm W x 900mm D
Product Code: FPD-10218



[LEC ULT531 Low Temperature 500 Litre Upright Freezer](#)

Capacity: 500 Litre
Config: Free Standing
Product Code: FPD-10219

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC CLT111 Low Temperature 113 Litre Chest Freezer



Product Code: FPD-10217

Dimensions: 1080mm H x 730mm W x 850mm D

Internal Dimensions: 450mm H x 450mm W x 560mm D

- -20°C to -40°C set temperature - Built for long term sample storage
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Special execution internal evaporator - To ensure temperature stability within the freezer with a +/- 1°C thermal gradient
- Spark-free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Colour: White
- Door Type: solid
- Door Style: hinged lid
- Lock
- Temperature Display on Panel
- Temperature Alarm
- Manual defrost
- Controller type: digital, external
- Refrigerant Type: R404a
- CFC Free
- UK 3 Pin Plug
- Voltage Supply: 230V
- Spark-free
- Aux Output
- Noise Level (dB): 52
- Energy Consumption kw/24h: 5

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ULT341 Low Temperature 322 Litre Upright Freezer



Product Code: FPD-10218

Dimensions: 1985mm H x 765mm W x 900mm D

Internal Dimensions: 1350mm H x 700mm W x 530mm D

- -20°C to -40°C set temperature - Ideal for medium term sample storage
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Special execution internal evaporator - To ensure temperature stability within the freezer with a +/- 1°C thermal gradient
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Colour: Silver/blue
- Door Type: solid
- Door Style: hinged
- Lock
- Temperature Display on Panel
- Temperature Alarm
- Number of Shelves: 2
- Auto or Manual defrost: manual
- Fan: noController type: digital, external
- Refrigerant Type: R404a
- CFC Free
- UK 3 Pin Plug
- Voltage Supply: 230V
- Spark-free
- Aux Output
- Noise Level (dB): 52
- Energy Consumption kw/24h: 6.5

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ULT531 Low Temperature 500 Litre Upright Freezer



Product Code: FPD-10219

Dimensions:

- -20°C to -40°C set temperature - Ideal for medium term sample storage
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Special execution internal evaporator - To ensure temperature stability within the freezer with a +/- 1°C thermal gradient
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- External Finish: blue, silver
- Door Type: solid
- Door Style: hinged
- Lock
- Temperature Display on Panel
- Temperature Alarm
- Number of Shelves: 2
- Manual defrost
- Controller type: digital, external
- Refrigerant Type: R404a
- CFC Free
- HFC Free
- UK 3 Pin Plug
- Voltage Supply: 230V
- Spark-free
- Aux Output
- Noise Level (dB): 52
- Energy Consumption kw/24h: 6.5

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Ultra Low Temperature Freezers

A variety of specialist Ultra Low Temperature Freezers Minus 40°C to Minus -85°C for the storage of chemicals, samples and other important matter.



[LEC CLT112 Low Temperature 113 Litre Chest Super Freezer](#)

Capacity: 113 Litre
Config: Chest
1080mm H x 730mm W x 850mm D
Product Code: FPD-10220



[LEC ULT342 Low Temperature 322 Litre Upright Super Freezer](#)

Capacity: 322 Litre
Config: Free Standing
1985mm H x 765mm W x 900mm D
Product Code: FPD-10221



[LEC ULT532 Low Temperature 500 Litre Upright Super Freezer](#)

Capacity: 500 Litre
Config: Free Standing
Product Code: FPD-10222

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC CLT112 Low Temperature 113 Litre Chest Super Freezer



Product Code: FPD-10220

Dimensions: 1080mm H x 730mm W x 850mm D

Internal Dimensions: 450mm H x 450mm W x 560mm D

- -40°C to -85°C set temperature - Built for long term sample storage
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Special execution internal evaporator - To ensure temperature stability within the freezer with a +/- 1°C thermal gradient
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Colour: White
- Door Type: solid
- Door Style: hinged lid
- Lock
- Temperature Display on Panel
- Temperature Alarm
- Manual defrost
- Controller type: digital, external
- Refrigerant Type: R404a/R508 mix
- CFC Free
- UK 3 Pin Plug
- Voltage Supply: 230V
- Spark-free
- Aux Output
- Noise Level (dB): 52
- Energy Consumption kw/24h: 12.4

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ULT342 Low Temperature 322 Litre Upright Super Freezer



Product Code: FPD-10221

Dimensions: 1985mm H x 765mm W x 900mm D

Internal Dimensions: 1350mm H x 700mm W x 530mm D

- -40°C to -85°C set temperature - Built for long term sample storage
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Special execution internal evaporator - To ensure temperature stability within the freezer with a +/- 1°C thermal gradient
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Colour: Silver/blue
- Door Type: solid
- Door Style: hinged
- Lock
- Temperature Display on Panel
- Temperature Alarm
- Number of Shelves: 2
- Manual defrost
- Controller type: digital, external
- Refrigerant Type: R404a/R508 mix
- CFC Free
- UK 3 Pin Plug
- Voltage Supply: 230V
- Spark-free
- Aux Output
- Noise Level (dB): 52
- Energy Consumption kw/24h: 16.2

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC ULT532 Low Temperature 500 Litre Upright Super Freezer



Product Code: FPD-10222

Dimensions:

- -40°C to -85°C set temperature - Built for long term sample storage
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the freezer
- Lock and two keys - for improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Special execution internal evaporator - to ensure temperature stability within the freezer with a +/- 1°C thermal gradient
- Spark free - There is nothing electrical inside the freezer so it is safe to store reactive substances in without the risk of explosion
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature starts to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- External Finish: blue, silver
- Colour: Silver/blue
- Door Type: solid
- Door Style: hinged
- Lock
- Temperature Display on Panel
- Temperature Alarm
- Number of Shelves: 2
- Manual defrost
- Controller type: digital, external
- Refrigerant Type: R404a/R508 mix
- CFC Free
- UK 3 Pin Plug
- Voltage Supply: 230V
- Spark-free
- Aux Output
- Noise Level (dB): 52
- Energy Consumption kw/24h: 21.5

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Laboratory Super Freezers

For laboratory and general use where storage at below normal freezing temperatures down to -40°C and/or where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



[Labcold RLBF1740 Superfreezer 50 Litres](#)

Capacity: 50 Litre
Config: Bench Top
580mm H x 550mm W x 620mm D

Product Code: FPD-04969



[Labcold RLBF1740/LK Superfreezer 50 Litres](#)

Capacity: 50 Litre
Config: Bench Top
580mm H x 550mm W x 620mm D

Product Code: FPD-04970



[Labcold RLHE0845 Superfreezer 237 Litres](#)

Capacity: 237 Litre
Config: Chest
860mm H x 1050mm W x 735mm D

Product Code: FPD-04972



[Labcold RLHE1145 Superfreezer 314 Litres](#)

Capacity: 314 Litre
Config: Chest
860mm H x 1305mm W x 735mm D

Product Code: FPD-04973



[Labcold RLHE1345 Superfreezer 375 Litres](#)

Capacity: 375 Litre
Config: Chest
860mm H x 1500mm W x 735mm D

Product Code: FPD-04974

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLBF1740/LK Superfreezer 50 Litres



Product Code: FPD-04970

Dimensions: 580mm H x 550mm W x 620mm D

- For laboratory and general use where storage at below normal freezing temperatures down to -40°C and/or where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity: 50 Litres (gross)
- Solid Door
- Door Lock
- Construction: White enamel exterior, moulded plastic interior
- Shelves supplied: 1, wire mesh
- Audible and visible high temperature alarm
- Digital temperature controller with display
- Sparkfree interior
- CFC and HCFC free
- Temperature set-point range: -30 to -40°C
- Operator controlled defrost
- Remote alarm connections available as optional extra

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLBF1740 Superfreezer 50 Litres



Product Code: FPD-04969

Dimensions: 580mm H x 550mm W x 620mm D

- For laboratory and general use where storage at below normal freezing temperatures down to -40°C and/or where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.
- Capacity: 50 Litres (gross)
- Solid Door
- No Door Lock
- Construction: White enamel exterior, moulded plastic interior
- Shelves supplied: 1, wire mesh
- Audible and visible high temperature alarm
- Digital temperature controller with display
- Sparkfree interior
- CFC and HCFC free
- Temperature set-point range: -30 to -40°C
- Operator controlled defrost
- Remote alarm connections available as optional extra

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLHE0845 Superfreezer 237 Litres



Product Code: FPD-04972

Dimensions: 860mm H x 1050mm W x 735mm D

- For laboratory and general use where there is a need to store at lower than normal freezing temperatures and/or there is concern that volatile, flammable materials or components might be stored now, or at some point in the future. Chest configuration provides improved temperature uniformity and stability with much reduced temperature rise following door openings.
- Capacity: 237 litres (gross)
- Construction: White enamel exterior, Aluminium interior, Additional insulation
- Baskets supplied: 1, wire mesh
- Digital display of temperature Sparkfree interior
- Door lock
- Easy-grip door handle
- CFC and HCFC free
- Temperature set-point range: -20 to -45°C
- Audible high temperature
- Operator controlled defrost

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RLHE1145 Superfreezer 314 Litres



Product Code: FPD-04973

Dimensions: 860mm H x 1305mm W x 735mm D

- For laboratory and general use where there is a need to store at lower than normal freezing temperatures and/or there is concern that volatile, flammable materials or components might be stored now, or at some point in the future. Chest configuration provides improved temperature uniformity and stability with much reduced temperature rise following door openings.
- Configuration: Chest
- Capacity: 314 Litres
- Sparkfree interior
- Temperature stability
- Lockable
- Alarm system
- CFC and HCFC free
- Digital temperature display
- Auto Defrost
- Two Shelves
- Operating temp: -45°C
- Door Type: Solid

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Compact Low Temperature Laboratory Freezers

Labcold Compact ULT freezers utilize a Japanese adaptation of a refrigeration technique tried and tested for many years in freezers operating at temperatures as low as -150°C. Using a single, commonly available, European compressor and fewer moving parts than more conventional ULT freezers these products will provide reliable storage at temperatures down to -83°C. This new and unique range is set to revolutionize the use of ultra low temperature freezers, particularly when -80°C storage is needed close to the work place.



[Labcold ULTF80 Compact Ultra Low Temperature Freezer 74 Litre](#)

Capacity: 74 Litre
Config: Chest
875mm H x 550mm W x 645mm D
Product Code: FPD-04990



[Labcold ULTF320 Compact Ultra Low Temperature Freezer 276 Litre](#)

Capacity: 276 Litre
Config: Chest
885mm H x 1260mm W x 695mm D
Product Code: FPD-04991



[Labcold ULTF420 Compact Ultra Low Temperature Freezer 383 Litre](#)

Capacity: 383 Litre
Config: Chest
885mm H x 1560mm W x 695mm D
Product Code: FPD-04992

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold ULTF80 Compact Ultra Low Temperature Freezer 74 Litre



Product Code: FPD-04990

Dimensions: 875mm H x 550mm W x 645mm D

- Labcold Compact ULT freezers utilize a Japanese adaptation of a refrigeration technique tried and tested for many years in freezers operating at temperatures as low as -150°C. Using a single, commonly available, European compressor and fewer moving parts than more conventional ULT freezers these products will provide reliable storage at temperatures down to -83°C. This new and unique range is set to revolutionize the use of ultra low temperature freezers, particularly when -80°C storage is needed close to the work place.
- With their compact external dimensions and quiet operation they are no more obtrusive than a standard lab refrigerator.
- Operating temperature -83°C
- 74 Litres
- Door lock
- Very quiet
- Low cost
- Unique application of proven technology
- Digital display with integral alarm
- Low heat output
- Voltage Supply: 230V/50Hz
- UK 3 Pin Plug

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold ULTF320 Compact Ultra Low Temperature Freezer 276 Litre

Product Code: FPD-04991

Dimensions: 885mm H x 1260mm W x 695mm D



- Labcold Compact ULT freezers utilize a Japanese adaptation of a refrigeration technique tried and tested for many years in freezers operating at temperatures as low as -150°C. Using a single, commonly available, European compressor and fewer moving parts than more conventional ULT freezers these products will provide reliable storage at temperatures down to -83°C. This new and unique range is set to revolutionize the use of ultra low temperature freezers, particularly when -80°C storage is needed close to the work place. With their compact external dimensions and quiet operation they are no more obtrusive than a standard lab refrigerator.
- Operating temperature -83°C
- Door lock
- Temperature Display
- Alarm
- Very quiet operation
- Low cost
- Unique application of proven technology
- Digital display with integral alarm

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold ULTF420 Compact Ultra Low Temperature Freezer 383 Litre

Product Code: FPD-04992

Dimensions: 885mm H x 1560mm W x 695mm D



- Labcold Compact ULT freezers utilize a Japanese adaptation of a refrigeration technique tried and tested for many years in freezers operating at temperatures as low as -150°C. Using a single, commonly available, European compressor and fewer moving parts than more conventional ULT freezers these products will provide reliable storage at temperatures down to -83°C. This new and unique range is set to revolutionize the use of ultra low temperature freezers, particularly when -80°C storage is needed close to the work place. With their compact external dimensions and quiet operation they are no more obtrusive than a standard lab refrigerator.
- Operating temperature -83°C
- Door lock
- Temperature Display
- Alarm
- Very quiet operation
- Low cost
- Unique application of proven technology
- Digital display with integral alarm

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

DuoCold Ultra Low Temperature Freezers

Designed to store materials safely at low temperatures, this NEW Labcold freezer features two compressors and is backed by a 5 year warranty comprising of three years parts and labour with maintenance visits with the option for a further 2 years parts only warranty, making this freezer a cost effective purchase. Featuring VIP technology, this freezer is designed to provide a large amount of storage whilst occupying as little laboratory space as possible.



[Labcold LULT0480 DuoCold Ultra Low Temperature Freezer 100 Litre](#)

Capacity: 100 Litre
Config: Chest
745mm H x 1095mm W x 750mm D

Product Code: FPD-10244



[Labcold LULT1480 DuoCold Ultra Low Temperature Freezer 388 Litre](#)

Capacity: 388 Litre
Config: Chest
2020mm H x 840mm W x 985mm D

Product Code: FPD-10245



[Labcold LULT2080 DuoCold Ultra Low Temperature Freezer 538 Litre](#)

Capacity: 538 Litre
Config: Chest
2020mm H x 840mm W x 985mm D

Product Code: FPD-10246

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold LULT0480 DuoCold Ultra Low Temperature Freezer 100 Litre

Product Code: FPD-10244

Dimensions: 745mm H x 1095mm W x 750mm D

Internal Dimensions: 515mm H x 455mm W x 450mm D



- The Labcold 100 litre ULT freezer features VIP technology and is designed to provide a large amount of storage whilst occupying as little laboratory space as possible.
- Featuring one energy saving Danfoss compressor this freezer can pull down to -86°C in an impressive 3-4 hours. Due to the accurate microprocessor controller, the set point can be adjusted by one degree increments making the freezer extremely versatile and adaptable.
- Fitted with a range of audible and visual alarms, and with a 72 hour controller battery back up, this freezer features all the safety features normally found on a much larger cabinet including password setting protection.
- For extra peace of mind, this freezer comes complete with a free 5 year warranty. Two free maintenance visits are included in the warranty making this an extremely
- Capacity: 100 Litre
- Configuration: Under bench
- Operating Temperature -86°C
- Interior Construction Stainless Steel
- Exterior Construction Powder Coated steel
- Door Style Solid, Right Hinged
- High Temperature Alarm, Low Temperature Alarm
- Mains Failure Alarm
- Low Battery Alarm
- Blocked Filter alarm
- Digital Temperature Display
- Stored Temperature Data (7 Days)
- Key Operated Door Lock
- Password Protected Temperature Settings
- Battery Back up (Display)
- Two Shelves
- Inner Door
- Vacuum Release Port
- Eco Friendly Refrigerant

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold LULT1480 DuoCold Ultra Low Temperature Freezer 388 Litre



Product Code: FPD-10245

Dimensions: 2020mm H x 840mm W x 985mm D

Internal Dimensions: 1200mm H x 500mm W x 600mm D

- The Labcold 388 litre ULT freezer is designed to provide secure storage of items at low temperature and is robustly constructed from steel with a stainless steel interior.
- Featuring two energy saving Danfoss compressors this freezer can pull down to -86° C in an impressive 4 hours. Due to the accurate microprocessor controller, the set point can be adjusted by one degree increments making the freezer extremely versatile and adaptable.
- Operating Temperature -86°C
- Interior Construction Stainless Steel
- Exterior Construction Powder Coated steel
- Door Style Solid, Right Hinged
- Open Door alarm
- High Temperature Alarm, Low Temperature Alarm
- Mains Failure Alarm
- Low Battery Alarm
- Blocked Filter alarm
- Digital Temperature Display
- Stored Temperature Data (7 Days)
- Key Operated Door Lock
- Password Protected Temperature Settings
- Battery Back up (Display)
- Three Shelves
- Four Inner Doors
- Four Castors
- Vacuum Release Port
- Eco Friendly Refrigerant

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold LULT2080 DuoCold Ultra Low Temperature Freezer 538 Litre



Product Code: FPD-10246

Dimensions: 2020mm H x 840mm W x 985mm D

Internal Dimensions: 1350mm H x 600mm W x 700mm D

- The Labcold 538 litre ULT freezer features VIP technology and is designed to provide a large amount of storage whilst occupying as little laboratory space as possible.
- Featuring two energy saving Danfoss compressors this freezer can pull down to -86° C in an impressive 4 hours. Due to the accurate microprocessor controller, the set point can be adjusted by one degree increments making the freezer extremely versatile and adaptable.
- Capacity (Litres) 538
- Configuration Freestanding
- Operating Temperature -86°C
- Interior Construction Stainless Steel
- Exterior Construction Powder Coated steel
- Door Style Solid
- Door Hinging Right
- Open Door alarm
- High Temperature Alarm, Low Temperature Alarm
- Mains Failure Alarm
- Low Battery Alarm
- Blocked Filter alarm
- Digital Temperature Display
- Stored Temperature Data (7 Days)
- Key Operated Door Lock
- Password Protected Temperature Settings
- Battery Back up (Display)
- Three Shelves
- Four Inner Doors
- Four Castors
- Vacuum Release Port
- Eco Friendly Refrigerant

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Ward Refrigeration



Ward Refrigerators

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Ward Fridge

This range provide the secure storage of foodstuffs within the ward environment. The newly improved range operate at +3°C with an available operating range of +1° to +5°C, all refrigerators come with lockable doors, external temperature controller and have a **two year warranty**.



[LEC WR207C Ward Refrigerator 82 Litre](#)

Capacity: 82 Litre
Config: Bench Top
660mm H x 502mm W x 540mm D
Product Code: FPD-09812



[LEC WSR151UK Ward Refrigerator 151 Litre](#)

Capacity: 151 Litre
Config: Under Bench
845mm H x 595mm W x 635mm D
Product Code: FPD-09813



[LEC WSR288UK Ward Refrigerator 288 Litre](#)

Capacity: 288 Litre
Config: Free Standing
1565mm H x 595mm W x 670mm D
Product Code: FPD-09814

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC WR207C Ward Refrigerator 82 Litre



Product Code: FPD-09812

Dimensions: 660mm H x 502mm W x 540mm D

Internal Dimensions: 540mm H x 419mm W x 315mm D

- +1°C to +5°C set temperature - Operates at +3°C to conform to the requirements of Food Hygiene Regulations
- External digital controller - For temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- Reversible door - For maximum versatility in sometimes limited working spaces
- Small countertop size - For limited space needing medical refrigeration
- External Finish: White painted steel
- Internal Finish: White plastic
- Type of plastic: HIPS / ABS Mix
- Colour: White/green
- Door Type: Solid - Hinged, reversible
- Door Open Depth: 976mm
- Temperature Display on Panel
- Two Wire Shelves
- Shelf Size: 440x275x6mm
- Automatic Defrost
- Climate Class: N
- Refrigerant Type: R134a
- CFC Free, HFC Free, HCFC Free
- UK 3 Pin Plug, Voltage Supply: 230V
- Noise Level (dB): <43
- Insulation Thickness: 30mm
- Insulation Type: Cyclopentane
- Energy Consumption kw/24h: 1.1

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC WSR151UK Ward Refrigerator 151 Litre



Product Code: FPD-09813

Dimensions: 845mm H x 595mm W x 635mm D

Internal Dimensions: 705mm H x 500mm W x 450mm D

- Set temperature - Operates at +3°C to keep food fresher for longer
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Roller castors - To make the fridge easier to move for cleaning and maintenance

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

LEC WSR288UK Ward Refrigerator 288 Litre



Product Code: FPD-09814

Dimensions: 1565mm H x 595mm W x 670mm D

Internal Dimensions: 1382mm H x 480mm W x 445mm D

- Set temperature - Operates at +3°C to keep food fresher for longer
- External digital controller accurate to 0.1°C - For precise temperature readings without having to open the door, so preventing unnecessary temperature fluctuations within the fridge
- Anti microbial handle - To help prevent the spread of bacteria and infection
- Lock and two keys - For improved safety, accountability and security over domestic refrigeration, a lockable fridge helps prevent unauthorised access to valuable and dangerous contents
- Fan assisted - To aid the flow of cool air around the fridge, eliminating temperature variations to keep the contents at a consistent temperature. Also speeds up temperature recovery after door opening
- Audible and visual high and low temperature alarms - To alert anybody nearby as soon as the temperature begins to go too high or low, so that action can be taken to rectify the problem and prevent damage to the contents
- R600a refrigerant - Environmentally friendly refrigerant with a reduced level of carbon dioxide given off, reduces running costs by reducing how hard the compressor has to work to keep the fridge cool
- Open door alarm - To prevent spoilage from raised internal temperature and illicit entry if the door is left open
- Reversible door - For maximum versatility in sometimes limited working spaces
- Interior Xw light - For better visibility of contents

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Portable Refrigeration & Vaccine Bags



Portable Refrigeration
Thermal Transportation Bags
Premium Vaccine Carriers

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Portable Refrigeration

A selection of products designed for transportation of Vaccines and Specimen Collection. The portable refrigerator comes with an in-car power cable to keep its contents at the desired temperature whilst in transit. The Thermal Transportation bags allow for transportation of its contents over short trips approx 3 hours And the premium Vaccine carriers can hold the temperature for up to Eight hours allowing for longer transit times.



[Labcold RPDF0012D Portable Medical Fridge](#)

£340.00

290mm H x 420mm W x 290mm D

Product Code: FPD-04387



[Shoreline SMP65 - 65 Litres Portable Medical Fridge](#)

£748.00

Capacity: 65 Litre
450mm H x 725mm W x 510mm D

Product Code: FPD-06632



[Shoreline SMP26 - 26 Litres Portable Medical Fridge](#)

£599.00

Capacity: 26 Litre
385mm H x 345mm W x 510mm D

Product Code: FPD-06254



[Shoreline SMP41 - 41 Litres Portable Medical Fridge](#)

£638.00

Capacity: 41 Litre
400mm H x 610mm W x 385mm D

Product Code: FPD-04389

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Labcold RPDF0012D Portable Medical Fridge



Product Code: FPD-04387

Dimensions: 290mm H x 420mm W x 290mm D

- Simple and versatile the Labcold Digital Portable Laboratory & Medical Transporter can both cool or heat and is designed for the transporting of samples, vaccines, pharmaceuticals or indeed anything that needs to be moved under temperature controlled conditions. Complete with an in-car power cable in addition to a mains cable, this transporter can keep contents at the desired temperature whether in transit or in situ.
- Configuration: Electronic Portable
- Capacity: 12 Litres
- Temperature range: -5 to +60°C
- Door Type: Top opening
- Two Compartments
- Weight: 6.2 Kg
- Supplied with mains and in car power cables
- Safe gas free peltier refrigeration
- Lightweight and easily transportable
- Digital set point display
- Easy digital temperature adjustment
- Tough plastic exterior
- **Please be aware this product is not suitable for the transport of blood or blood products and is not designed to be run for extended periods of time.**

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SMP26 - 26 Litres Portable Medical Fridge



Product Code: FPD-06254

Dimensions: 385mm H x 345mm W x 510mm D

Internal Dimensions: 255mm H x 275mm W x 435mm D

- 12 Volt DC (fitted with Car cigarette lighter connection for mobile operation)
- or 220/110 Volt AC (fitted with 3 pin plug for mains power operation)
-
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- CFC & Ammonia free
- Finished in durable-moulded exterior
- 16 Kgs
- 33 Watts
- Max. Ambient +32°C

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP™

Shoreline SMP41 - 41 Litres Portable Medical Fridge



Product Code: FPD-04389

Dimensions: 400mm H x 610mm W x 385mm D

Internal Dimensions: 260mm H x 525mm W x 290mm D

- For the safe storage of Temperature Sensitive Products during transit.
- Capacity: 41 litres
- 12 Volt DC
- Maintains internal air temperature between +2°C and +8°C. External digital temperature display with MAX/MIN memory for continuous monitoring. Forced-air cooling for temperature stability and rapid temperature recovery after lid openings. Audio/visual alarm signal on temperature deviation.
- Finished in durable-moulded exterior
- 22 Kgs.

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Shoreline SMP65 - 65 Litres Portable Medical Fridge



Product Code: FPD-06632

Dimensions: 450mm H x 725mm W x 510mm D

Internal Dimensions: 325mm H x 610mm W x 360mm D

- 12 Volt DC (fitted with Car cigarette lighter connection for mobile operation)
- or 220/110 Volt AC (fitted with 3 pin plug for mains power operation)
- 65 litres Capacity
- +2°C to +8°C
- Forced-Air Cooling
- Digital Temperature Display
- Audio / Visual Temperature Alarm
- MAX/MIN Temperature Memory
- CFC & Ammonia free
- Finished in durable-moulded exterior
- 24 Kgs
- 42 Watts
- Max. Ambient +32°C
- Prices are Exclusive of VAT
- Prices Includes Delivery Mainland UK

Delivery

Typical Lead Time: 10-14 Days from receipt of payment.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP

Thermal Transportation Bags & Premium Vaccine Carriers

Our standard range of 10, 20, and 30 litre carry bags are ideal for vaccination programs lasting up to 8 hours. Our bags lose 0.74°C per hour without the use of coolants giving you 8.1 hours if the samples were put in at 2°C and the lid unopened! We do however advise the use of 5°C cool packs to help maintain temperature. WE advise a minimum of a cool pack on the top and bottom of the samples but ideally one on each side also. If you are opening and closing the lid on the bag repeatedly, use a layering system in the bag by incorporating our thermal separators in between each layer of vaccines. This reduces the impact of the ambient air getting in to the bag on opening. Also the lid should be opened the bare minimum to extract the vaccines and for the least amount of time!

The 3 Inch and 2 inch thermal vaccine carriers have been designed to maintain a temperature of between a 2° to 8°C for up to 8 hours, especially for use with vaccination programmes or specimen collection. The inner lid reduces thermal leakage helping to maintain temperature parameters despite repeated opening of the bag. Developed using NASA technology Polartherm™ is made up of an inner reflective layer, a hollow fibre filler and a waterproof woven PVC outer, the result of which is a material that is both tough and lightweight with excellent thermal control properties.



[Two Inch Vaccine Thermal Carrier](#)

£298.00

Config: 2 Inch Deep Internal Chamber
240mm H x 320mm W x 280mm D

Product Code: FPD-10140



[Three Inch Vaccine Thermal Carrier](#)

£298.00

Config: 3 Inch Deep Internal Chamber
240mm H x 320mm W x 280mm D

Product Code: FPD-10141



[10 Litre Thermal Carry Bag with 1 Separator](#)

£161.59

Capacity: 10 Litre
220mm H x 400mm W x 280mm D

Product Code: FPD-06703



[20 Litre Thermal Carry Bag with 2 Separators](#)

£168.00

Capacity: 20 Litre
330mm H x 400mm W x 280mm D

Product Code: FPD-06704



[30 Litre Thermal Carry Bag with 3 Separators](#)

£201.96

Capacity: 30 Litre
360mm H x 400mm W x 340mm D

Product Code: FPD-06705



[30 Litre Thermal Carry Bag with Trolley and 3 Separators](#)

£241.56

Capacity: 30 Litre
360mm H x 400mm W x 340mm D

Product Code: FPD-06706



[A4 Thermal Vaccine Wallet](#)

£82.28

260mm H x 360mm W

Product Code: FPD-06707

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Two Inch Vaccine Thermal Carrier



Product Code: FPD-10140

Dimensions: 240mm H x 320mm W x 280mm D

Internal Dimensions: 50mm H x 200mm W x 180mm D

- The 3 Inch and 2 inch thermal vaccine carriers have been designed to maintain a temperature of between a 2° to 8°C for up to 8 hours, especially for use with vaccination programmes or specimen collection. The inner lid reduces thermal leakage helping to maintain temperature parameters despite repeated opening of the bag. Developed using NASA technology Polartherm? is made up of an inner reflective layer, a hollow fibre filler and a waterproof woven PVC outer, the result of which is a material that is both tough and lightweight with excellent thermal control properties. The system comes with 4 thermal bottles which should be pre conditioned before use as follows:
- 2 x StrurdeeSeal?2 Bottles (Black Cap) - Pre condition at 5°C for a minimum of 48 Hours prior to assembly
- 2 x StrurdeeSeal?1 Bottles (White Cap) - Precondition at -20°C for a minimum of 48 Hours prior to assembly
- Any clinical samples should be preconditioned at 5°C for 8 hours prior to assembly
- The -20°C Bottles should be placed directly against the lid and the bottom of the carrier and the 5°C should be placed against the -20°C bottle in the lid and at the bottom of the carrier.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Three Inch Vaccine Thermal Carrier



Product Code: FPD-10141

Dimensions: 240mm H x 320mm W x 280mm D

Internal Dimensions: 75mm H x 200mm W x 180mm D

- The 3 Inch and 2 inch thermal vaccine carriers have been designed to maintain a temperature of between a 2° to 8°C for up to 8 hours, especially for use with vaccination programmes or specimen collection. The inner lid reduces thermal leakage helping to maintain temperature parameters despite repeated opening of the bag. Developed using NASA technology Polartherm? is made up of an inner reflective layer, a hollow fibre filler and a waterproof woven PVC outer, the result of which is a material that is both tough and lightweight with excellent thermal control properties. The system comes with 4 thermal bottles which should be pre conditioned before use as follows:
- 2 x StrurdeeSeal?2 Bottles (Black Cap) - Pre condition at 5°C for a minimum of 48 Hours prior to assembly
- 2 x StrurdeeSeal?1 Bottles (White Cap) - Precondition at -20°C for a minimum of 48 Hours prior to assembly
- Any clinical samples should be preconditioned at 5°C for 8 hours prior to assembly
- The -20°C Bottles should be placed directly against the lid and the bottom of the carrier and the 5°C should be placed against the -20°C bottle in the lid and at the bottom of the carrier.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

10 Litre Thermal Carry Bag with 1 Separator



Product Code: FPD-06703

Dimensions: 220mm H x 400mm W x 280mm D

- Weight empty: 0.4648 kgs
- External surface material: Cotton coated PVC
- Internal lining material: Cotton coated PVC
- Insulation material: Patented Polartherm material
- Insulation thickness: 45mm
- Internal dimensions: width: 350mm height: 170mm depth: 230mm
- Lid type and fixings: Zip
- Number of icepacks supplied: None
- Cold life without opening: Over 24hrs, temperature increase at an approximate rate of 0.73°C/hr from 3.5°C at an ambient temperature of 30°C. Please ask for copy of independent tests.
- Thermal Conductance (W/m².K): 1.01 ± 0.1
- Thermal Conductivity (W/m.K): 0.05 ± 0.005
- * Recommended use with 2 x One Litre Gel Packs

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

20 Litre Thermal Carry Bag with 2 Separators



Product Code: FPD-06704

Dimensions: 330mm H x 400mm W x 280mm D

- Weight empty: 0.9297 kgs
- External surface material: Cotton coated PVC
- Internal lining material: Cotton coated PVC
- Insulation material: Patented Polartherm material
- Insulation thickness: 45mm
- Internal dimensions: width: 350mm height: 280mm depth: 230mm
- Lid type and fixings: Zip
- Number of icepacks supplied: None
- Cold life without opening: Over 24hrs, temperature increase at an approximate rate of 0.73°C/hr from 3.5°C at an ambient temperature of 30°C. Please ask for copy of independent tests.
- Thermal Conductance (W/m².K): 1.01 ± 0.1
- Thermal Conductivity (W/m.K): 0.05 ± 0.005
- * Recommended use with 4 x One Litre Gel Packs

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

30 Litre Thermal Carry Bag with Trolley and 3 Separators



Product Code: FPD-06706

Dimensions: 360mm H x 400mm W x 340mm D

- Weight empty: 1.3945 kgs
- External surface material: Cotton coated PVC
- Internal lining material: Cotton coated PVC
- Insulation material: Patented Polartherm material
- Insulation thickness: 45mm
- External dimensions: width: 420mm height: 360mm depth: 340mm
- Lid type and fixings: Zip
- Number of icepacks supplied: None
- Cold life without opening: Over 24hrs, temperature increase at an approximate rate of 0.73°C/hr from 3.5°C at an ambient temperature of 30°C. Please ask for copy of independent tests.
- Thermal Conductance (W/m².K): 1.01 ± 0.1
- Thermal Conductivity (W/m.K): 0.05 ± 0.005
- * Recommended use with 6 x One Litre Gel Packs

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

30 Litre Thermal Carry Bag with 3 Separators



Product Code: FPD-06705

Dimensions: 360mm H x 400mm W x 340mm D

- Weight empty: 1.3945 kgs
- External surface material: Cotton coated PVC
- Internal lining material: Cotton coated PVC
- Insulation material: Patented Polartherm material
- Insulation thickness: 45mm
- Internal dimensions: width: 380mm height: 320mm depth: 300mm
- Lid type and fixings: Zip
- Number of icepacks supplied: None
- Cold life without opening: Temperature increase at an approximate rate of 0.67°C/hr when the internal/external temperature difference is over 20°C. (Minimum) Please ask for copy of independent tests.
- Thermal Conductance (W/m².K): 1.01 ± 0.1
- Thermal Conductivity (W/m.K): 0.05 ± 0.005
-
- * Recommended use with 6 x One Litre Gel Packs

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

A4 Thermal Vaccine Wallet



Product Code: FPD-06707

Dimensions: 260mm H x 360mm W

- Vaccine storage capacity: 2 ltrs
- Weight empty: 0.300 kgs
- External surface material: Cotton coated PVC
- Internal lining material: Cotton coated PVC
- Insulation material: Patented Polartherm material
- Insulation thickness: 45mm
- Internal dimensions: width: 310mm height: 210mm
- Vaccine storage dimensions: width: 310mm height: 210mm (max)
- Lid type and fixings: Velcro
- Number of icepacks supplied: None
- Cold life without opening: 500ml of ice enclosed in the pouch will remain at 0°C for 18 hours in a 30°C ambient.
- Thermal Conductance (W/m².K): 1.01 ± 0.1
- Thermal Conductivity (W/m.K): 0.05 ± 0.005
- * Recommended use with 1 x Half Litre Gel Packs

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Laboratory Incubators



Laboratory Cooled Incubators
Laboratory Refrigerated Cooled Incubators

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Laboratory Incubators

Cooled incubators are ideal for microbiology or any application where precise and consistent temperatures at or below room temperature are demanded. Available in a range of sizes from a small benchtop model through to a large 600 litre unit constructed from mild steel with stainless steel interior, these new cooled incubators from Labcold all feature fan air circulation for even temperature distribution and microprocessor controlled temperature regulation for improved accuracy and stability.



[RLCG0100 Labcold Cooled Incubator with Glass Door](#)

Capacity: 88 Packs
Config: Bench Top/Wall Mounted
810mm H x 480mm W x 530mm D

Product Code: FPD-10110



[RLCG01502 Labcold Cooled Incubator with Glass Door](#)

Capacity: 150 Packs
Config: Under Bench
820mm H x 600mm W x 600mm D

Product Code: FPD-10111



[RLSD01502 Labcold Cooled Incubator with Solid Door](#)

Capacity: 150 Packs
Config: Under Bench
820mm H x 600mm W x 600mm D

Product Code: FPD-10112



[RLSD0300 Labcold Refrigerated Cooled Incubator - Solid Door](#)

Capacity: 300 Packs
Config: Free Standing
1487mm H x 600mm W x 650mm D

Product Code: FPD-10113



[RLCG0600 Labcold Refrigerated Cooled Incubator - Glass Door](#)

Capacity: 600 Packs
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-10114



[RLSD0600 Labcold Refrigerated Cooled Incubator - Solid Door](#)

Capacity: 600 Packs
Config: Free Standing
1980mm H x 720mm W x 820mm D

Product Code: FPD-10115

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

RLCG0100 Labcold Cooled Incubator with Glass Door



Product Code: FPD-10110

Dimensions: 810mm H x 480mm W x 530mm D

- Benchtop controlled temperature incubator for maintaining samples at temperatures accurately controlled within the range +5 to +40°C.
- Capacity: 88 litres (gross)
- Door Type: Glass
- Construction: White enamel exterior, Stainless steel interior, Double glazed door.
- Shelves supplied: 2; Plastic coated wire mesh
- Microprocessor temperature control with digital display
- Temperature precision: better than $\pm 0.5^{\circ}\text{C}$ at +5°C
- Fan circulated interior
- CFC and HFC free
- Manual defrost with condensate drain
- Maximum ambient temperature: +32°C

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

RLSD01502 Labcold Cooled Incubator with Solid Door



Product Code: FPD-10112

Dimensions: 820mm H x 600mm W x 600mm D

- Temperature controlled incubator for maintaining samples at temperatures accurately within the range +5 to +40°C.
- Capacity: 150 litres (gross)
- Construction: White enamel exterior, moulded plastic interior.
- Shelves supplied: 3; Plastic coated wire mesh
- Microprocessor temperature control with digital display
- Temperature set-point range: +5 to +40°C
- Temperature precision: better than $\pm 0.5^{\circ}\text{C}$ at +5°C
- Fan circulated interior
- Integral Handle
- Door Lock
- CFC and HFC free
- Manual defrost with condensate drain
- Maximum ambient temperature: +32°C

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

RLCG01502 Labcold Cooled Incubator with Glass Door



Product Code: FPD-10111

Dimensions: 820mm H x 600mm W x 600mm D

- Temperature controlled incubator for maintaining samples at temperatures accurately within the range +5 to +40°C.
- Capacity: 150 litres (gross)
- Construction: White enamel exterior, Moulded plastic interior, Double glazed door.
- Shelves supplied: 3; Plastic coated wire mesh
- Microprocessor temperature control with digital display
- Temperature set-point range: +5 to +40°C
- Temperature precision: better than $\pm 0.5^{\circ}\text{C}$ at +5°C
- Fan circulated interior
- Integral Handle
- Door Lock
- CFC and HFC free
- Manual defrost with condensate drain
- Maximum ambient temperature: +32°C

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

RLSD0300 Labcold Refrigerated Cooled Incubator - Solid Door



Product Code: FPD-10113

Dimensions: 1487mm H x 600mm W x 650mm D

- This 300 litre cooled incubator is ideal for microbiology or any application where precise and consistent temperatures at or below room temperature are required. With a robust painted steel interior and exterior and fully adjustable shelves, this incubator is versatile and built to last.
- Capacity (Litres) 300
- Config: Freestanding
- Weight 75kg
- Energy Consumption 0.9kWh/24hr
- Operating Temperature +5°C to +40°C
- Interior Construction Painted Steel
- Exterior Construction Powder Coated steel
- Door Style Solid
- Door Hinging Right
- Lock Manual Key
- Shelves 3
- Shelf dimension (WxD) 495 x 460mm
- Spare shelves available
- Temperature display
- Microprocessor controlled temperature regulation
- Fan assisted temperature control
- Access Port
- Eco Friendly Refrigerant

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

RLSD0600 Labcold Refrigerated Cooled Incubator - Solid Door



Product Code: FPD-10115

Dimensions: 1980mm H x 720mm W x 820mm D

- Temperature controlled incubator for maintaining samples at temperatures accurately within the range +5 to +40°C.
- Capacity: 600 litres (gross)
- Construction: White and blue enamel exterior, stainless steel interior.
- Cabinet mounted on castors.
- Shelves supplied: 4; Stainless steel
- Microprocessor temperature control with digital display
- Temperature set-point range: +5 to +40°C
- Temperature precision: better than $\pm 0.5^{\circ}\text{C}$ at +5°C
- Fan circulated interior
- Inner clear view doors
- Integral Handle
- CFC and HFC free
- Manual defrost
- Maximum ambient temperature: +32°C

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

RLCG0600 Labcold Refrigerated Cooled Incubator - Glass Door



Product Code: FPD-10114

Dimensions: 1980mm H x 720mm W x 820mm D

- Temperature controlled incubator for maintaining samples at temperatures accurately within the range +5 to +40°C.
- Capacity: 600 litres (gross)
- Construction: White and blue enamel exterior, stainless steel interior.
- Cabinet mounted on castors.
- Shelves supplied: 4; Stainless steel
- Double glazed door
- Microprocessor temperature control with digital display
- Temperature set-point range: +5 to +40°C
- Temperature precision: better than $\pm 0.5^{\circ}\text{C}$ at +5°C
- Fan circulated interior
- Inner clear view doors
- Integral Handle
- CFC and HFC free
- Manual defrost
- Maximum ambient temperature: +32°C

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Bloodbanks



Solid Door Bloodbanks
Glass Door Bloodbanks

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Bloodbank Refrigeration

Lec Medical bloodbanks have been designed in conjunction with bloodbank centres and we are proud of the many features they provide. Each bloodbank is built to order in the UK and can be delivered directly to your haematology department. All bloodbanks come with a 48 hour battery back-up for the controller, an energy efficient LED light and auto defrost as standard.



[Blood Bank 48 Pack with Glass Door](#)

Capacity: 78 Packs
780mm H x 522mm W x 600mm D
Product Code: FPD-09661



[Blood Bank 48 Pack with Solid Door](#)

Capacity: 78 Packs
780mm H x 522mm W x 600mm D
Product Code: FPD-09662



[Blood Bank 48 Pack with Glass Door](#)

Capacity: 78 Packs
830mm H x 520mm W x 600mm D
Product Code: FPD-09663



[Blood Bank 48 Pack with Solid Door](#)

Capacity: 78 Packs
830mm H x 520mm W x 600mm D
Product Code: FPD-09664



[Blood Bank 72 Pack with Glass Door](#)

Capacity: 110 Packs
1200mm H x 520mm W x 605mm D
Product Code: FPD-09665



[Blood Bank 72 Pack with Solid Door](#)

Capacity: 110 Packs
1200mm H x 520mm W x 605mm D
Product Code: FPD-09666



[Blood Bank 360 Pack with Glass Door](#)

Capacity: 350 Packs
1935mm H x 680mm W x 865mm D
Product Code: FPD-09667



[Blood Bank 360 Pack with Solid Door](#)

Capacity: 350 Packs
1935mm H x 680mm W x 865mm D
Product Code: FPD-09668

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Bloodbank Refrigeration

Lec Medical bloodbanks have been designed in conjunction with bloodbank centres and we are proud of the many features they provide. Each bloodbank is built to order in the UK and can be delivered directly to your haematology department. All bloodbanks come with a 48 hour battery back-up for the controller, an energy efficient LED light and auto defrost as standard.



[Blood Bank 432 Pack with Glass Door](#)

Capacity: 439 Packs
1980mm H x 680mm W x 850mm D
Product Code: FPD-09669



[Blood Bank 432 Pack with Solid Door](#)

Capacity: 439 Packs
1980mm H x 680mm W x 850mm D
Product Code: FPD-09670



[Blood Bank 864 Pack with Glass Door](#)

Capacity: 960 Packs
1980mm H x 1340mm W x 850mm D
Product Code: FPD-09671



[Blood Bank 864 Pack with Solid Door](#)

Capacity: 960 Packs
1980mm H x 1340mm W x 850mm D
Product Code: FPD-09672

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 48 Pack with Solid Door



Product Code: FPD-09662

Dimensions: 780mm H x 522mm W x 600mm D

Internal Dimensions: 352mm H x 420mm W x 440mm D

- Single Solid Door
- Configuration: Countertop
- Refrigerant: R600a
- Stores 48 Blood Bags (250ml fill) or 28 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max./Min. Ambient Working Temperature: 40°C/10°C
- Digital lock with 2 key fobs
- Two Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish: Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Castors: None, One adjustable foot, Three fixed feet
- Tangential Fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant), Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN
- Energy Consumption 1.4 kWh/24h

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 48 Pack with Glass Door



Product Code: FPD-09661

Dimensions: 780mm H x 522mm W x 600mm D

Internal Dimensions: 352mm H x 420mm W x 440mm D

- Single Solid Door
- Configuration: Countertop
- Refrigerant: R600a
- Stores 48 Blood Bags (250ml fill) Or 28 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max/Min Ambient Working Temp: 40°C/10°C
- Digital lock with 2 key fobs
- Two Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Castors: None, One adjustable foot, Three fixed feet
- Tangential Fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant)
- Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN
- Energy Consumption 1.4 kWh/24h

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 72 Pack with Solid Door



Product Code: FPD-09666

Dimensions: 1200mm H x 520mm W x 605mm D

Internal Dimensions: 602mm H x 440mm W x 440mm D

- Single Solid Door
- Configuration: Freestanding
- Refrigerant: R600a
- Stores 72 Blood Bags (250ml fill) or 42 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max./Min. Ambient Working Temperature: 43°C/10°C
- Digital lock with 2 key fobs
- Three Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish: Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Four castors (front two lockable), Two adjustable feet
- Tangential Fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant), Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN
- Energy Consumption 1.5 kWh/24h

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 72 Pack with Glass Door



Product Code: FPD-09665

Dimensions: 1200mm H x 520mm W x 605mm D

Internal Dimensions: 602mm H x 440mm W x 440mm D

- Single Glass Door
- Configuration: Freestanding
- Refrigerant: R600a
- Stores 72 Blood Bags (250ml fill) or 42 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max./Min. Ambient Working Temperature: 40°C/10°C
- Digital lock with 2 key fobs
- Three Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish: Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Four castors (Front Two castors lockable), Two Adjustable feet
- Tangential fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant), Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN
- Energy Consumption 1.4 kWh/24h

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 360 Pack with Solid Door



Product Code: FPD-09668

Dimensions: 1935mm H x 680mm W x 865mm D

Internal Dimensions: 1000mm H x 550mm W x 640mm D

- Single Solid Door
- Configuration: Freestanding
- Refrigerant: R600a
- Stores 360 Blood Bags (250ml fill) or 210 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max./Min. Ambient Working Temperature: 43°C/10°C
- Digital lock with 2 key fobs
- Five Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish: Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Four castors (front two lockable), Two adjustable feet
- Axial Fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant), Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN
- Energy Consumption 3.4 kWh/24h

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 360 Pack with Glass Door



Product Code: FPD-09667

Dimensions: 1935mm H x 680mm W x 865mm D

Internal Dimensions: 1000mm H x 550mm W x 640mm D

- Single Glass Door
- Configuration: Freestanding
- Refrigerant: R600a
- Stores 360 Blood Bags (250ml fill) or 210 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max./Min. Ambient Working Temperature: 40°C/10°C
- Digital lock with 2 key fobs
- Five Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish: Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Four castors (front two lockable), Two adjustable feet
- Axial Fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant), Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN
- Energy Consumption 3.4 kWh/24h

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 432 Pack with Solid Door



Product Code: FPD-09670

Dimensions: 1980mm H x 680mm W x 850mm D

Internal Dimensions: 1225mm H x 560mm W x 640mm D

- Single Solid Door
- Configuration: Freestanding
- Refrigerant: R134a
- Stores 432 Blood Bags (250ml fill) or 252 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max./Min. Ambient Working Temperature: TBC/10°C
- Digital lock with 2 key fobs
- Dual refrigeration
- Six Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish: Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Four castors (front two lockable)
- Two adjustable feet
- Axial Fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant), Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 432 Pack with Glass Door



Product Code: FPD-09669

Dimensions: 1980mm H x 680mm W x 850mm D

Internal Dimensions: 1225mm H x 560mm W x 640mm D

- Single Glass Door
- Configuration: Freestanding
- Refrigerant: R134a
- Stores 432 Blood Bags (250ml fill) or 252 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max./Min. Ambient Working Temperature: TBC/10°C
- Digital lock with 2 key fobs
- Dual refrigeration
- Six Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish: Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Four castors (front two lockable), Two adjustable feet
- Axial Fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant), Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 864 Pack with Solid Door



Product Code: FPD-09672

Dimensions: 1980mm H x 1340mm W x 850mm D

Internal Dimensions: 1225mm H x 1225mm W x 640mm D

- Double Solid Door
- Configuration: Freestanding
- Refrigerant: R134a
- Stores 864 Blood Bags (250ml fill) or 505 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max./Min. Ambient Working Temperature: TBC/10°C
- Digital lock with 2 key fobs
- Dual refrigeration
- 12 Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish: Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Four castors (front two lockable), Two adjustable feet
- Axial Fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant), Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Blood Bank 864 Pack with Glass Door



Product Code: FPD-09671

Dimensions: 1980mm H x 1340mm W x 850mm D

Internal Dimensions: 1225mm H x 1225mm W x 640mm D

- Double Glass Door
- Configuration: Freestanding
- Refrigerant: R134a
- Stores 864 Blood Bags (250ml fill) or 505 Blood Bags (450ml fill)
- Working Temperature: 4°C, Load Temperature (bag temp): 2°C to 6°C
- Max./Min. Ambient Working Temperature: TBC/10°C
- Digital lock with 2 key fobs
- Dual refrigeration
- 12 Drawers
- Drawer construction: Stainless steel with clear perspex front
- External construction: Antic bacteria White powder coated, Zintec steel
- Internal finish: Stainless steel (grade 304)
- Internal Light (LED)
- Core temp alarm (Audible/Visual)
- Air temp alarm (Audible/Visual)
- Door open alarm (Audible/Visual)
- Mains failure alarm (Audible/Visual)
- Back up battery low alarm (Visual)
- Compressor failure alarm
- Load, Air and defrost sensor fail alarm
- Digital Display (LCD)
- Controller battery Back-Up
- Medical devices Certification
- BS4376 Certification
- RS485: None
- USB
- Chart recorder
- Four castors (front two lockable), Two adjustable feet
- Axial Fan
- Automatic defrost
- CFC free (Foam), HFC free (Refrigerant), HCFC free (Refrigerant), Amonia free
- Voltage: 230v 50Hz/60Hz, Plug Type: UK 3 PIN

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Icemakers & Water Dispensers



Ice/Water Dispensers
Flake/Nugget Ice Machines

FPDGROUP Ltd

www.pharmacy-fridge.co.uk

Icemakers & Water Dispensers

Compact ice/water dispenser - Perfect for use in busy areas where people need access to an ongoing fresh water and ice supply, such as a hospital ward. Hygienic ice production - Ice moves straight to an airtight storage hopper as soon as it is produced, so it has no contact with light or air until it is dispensed, keeping its hygienic integrity.



[Hoshizaki Countertop Ice/Water Dispenser \(60kg/24h\)](#)

Capacity: 60 kg/24h
Config: Free Standing
720mm H x 350mm W x 480mm D

Product Code: FPD-10230



[Hoshizaki Freestanding Ice/Water Dispenser \(32kg/24h\)](#)

Capacity: 32 kg/24h
Config: Free Standing
1600mm H x 350mm W x 500mm D

Product Code: FPD-10232



[Hoshizaki Countertop Ice/Water Dispenser \(125kg/24h\)](#)

Capacity: 152 kg/24h
Config: Free Standing
855mm H x 350mm W x 585mm D

Product Code: FPD-10231

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Icemakers & Water Dispensers

Compact ice/water dispenser - Perfect for use in busy areas where people need access to an ongoing fresh water and ice supply, such as a hospital ward. Hygienic ice production - Ice moves straight to an airtight storage hopper as soon as it is produced, so it has no contact with light or air until it is dispensed, keeping its hygienic integrity.



[Hoshizaki Freestanding Flake/Nugget Ice Machine \(85kg/24h\)](#)

Capacity: 85 kg/24h
Config: Free Standing
800mm H x 640mm W x 600mm D

Product Code: FPD-10233



[Hoshizaki Freestanding Flake/Nugget Ice Machine \(65kg/24h\)](#)

Capacity: 65 kg/24h
Config: Free Standing
800mm H x 640mm W x 600mm D

Product Code: FPD-10234



[Hoshizaki Freestanding Flake/Nugget Ice Machine \(125kg/24h\)](#)

Capacity: 125 kg/24h
Config: Free Standing
800mm H x 640mm W x 600mm D

Product Code: FPD-10235



[Hoshizaki Freestanding Flake/Nugget Ice Machine \(110kg/24h\)](#)

Capacity: 110 kg/24h
Config: Free Standing
800mm H x 640mm W x 600mm D

Product Code: FPD-10236



[Hoshizaki Freestanding Flake/Nugget Ice Machine \(125kg/24h\)](#)

Capacity: 125 kg/24h
Config: Free Standing
800mm H x 940mm W x 600mm D

Product Code: FPD-10237



[Hoshizaki Freestanding Flake/Nugget Ice Machine \(110kg/24h\)](#)

Capacity: 110 kg/24h
Config: Free Standing
800mm H x 940mm W x 600mm D

Product Code: FPD-10238

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Hoshizaki Countertop Ice/Water Dispenser (60kg/24h)



Product Code: FPD-10230

Dimensions: 720mm H x 350mm W x 480mm D

- Maintenance contract available - Consists of two visits per year at six month intervals to clean, sanitise and check operation, to help prevent breakdown and cut down on staff cleaning times
- 2 years parts and labour warranty - To reassure you of the high quality of our products, and give you added peace of mind should anything happen with these important products
- Compact ice/water dispenser - Perfect for use in busy areas where people need access to an ongoing fresh water and ice supply, such as a hospital ward
- Hygienic ice production - Ice moves straight to an airtight storage hopper as soon as it is produced, so it has no contact with light or air until it is dispensed, keeping its hygienic integrity
- Colour: Grey
- Climate Class: +5 to +40°C
- Refrigerant Type: R134a
- UK 3 Pin Plug
- Voltage Supply: 220-240V
- Energy Consumption kw/24h: 6.72
- Freezing Capacity KG per 24H: 60

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Hoshizaki Countertop Ice/Water Dispenser (125kg/24h)

Product Code: FPD-10231

Dimensions: 855mm H x 350mm W x 585mm D



- Maintenance contract available - Consists of two visits per year at six month intervals to clean, sanitise and check operation, to help prevent breakdown and cut down on staff cleaning times
- 2 years parts and labour warranty - To reassure you of the high quality of our products, and give you added peace of mind should anything happen with these important products
- Compact ice/water dispenser - Perfect for use in busy areas where people need access to an ongoing fresh water and ice supply, such as a hospital ward
- Hygienic ice production - Ice moves straight to an airtight storage hopper as soon as it is produced, so it has no contact with light or air until it is dispensed, keeping its hygienic integrity
- Colour: grey
- Climate Class: +5 to +40°C
- Refrigerant Type: R134
- UK 3 Pin Plug
- Voltage Supply: 220-240V
- Energy Consumption kw/24h: 14.5
- Freezing Capacity KG per 24H: 125

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Hoshizaki Freestanding Ice/Water Dispenser (32kg/24h)



Product Code: FPD-10232

Dimensions: 1600mm H x 350mm W x 500mm D

- Maintenance contract available - Consists of two visits per year at six month intervals to clean, sanitise and check operation, to help prevent breakdown and cut down on staff cleaning times
- 2 years parts and labour warranty - To reassure you of the high quality of our products, and give you added peace of mind should anything happen with these important products
- Compact ice/water dispenser - Perfect for use in busy areas where people need access to an ongoing fresh water and ice supply, such as a hospital ward
- Hygienic ice production - Ice moves straight to an airtight storage hopper as soon as it is produced, so it has no contact with light or air until it is dispensed, keeping its hygienic integrity
- Colour: Grey
- Climate Class: +1 to +40°C
- Refrigerant Type: R134a
- UK 3 Pin Plug
- Voltage Supply: 220-240V
- Energy Consumption kw/24h: 7.68
- Freezing Capacity KG per 24H: 32

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Hoshizaki Freestanding Flake/Nugget Ice Machine (85kg/24h)



Product Code: FPD-10233

Dimensions: 800mm H x 640mm W x 600mm D

- Maintenance contract available - Consists of two visits per year at six month intervals to clean, sanitise and check operation, to help prevent breakdown and cut down on staff cleaning times
- 2 years parts and labour warranty - To reassure you of the high quality of our products, and give you added peace of mind should anything happen with these important products
- Choice of flake or nugget models - To suit individual requirements
- Ice is compressed during production - Gives longer lasting cooling for transportation of blood and organs
- Colour: grey
- Climate Class: +5 to +40°C
- Refrigerant Type: R134a
- UK 3 Pin Plug
- Voltage Supply: 230V
- Energy Consumption kw/24h: 7.2
- Freezing Capacity KG per 24H: 85

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Hoshizaki Freestanding Flake/Nugget Ice Machine (65kg/24h)



Product Code: FPD-10234

Dimensions: 800mm H x 640mm W x 600mm D

- Maintenance contract available - Consists of two visits per year at six month intervals to clean, sanitise and check operation, to help prevent breakdown and cut down on staff cleaning times
- 2 years parts and labour warranty - To reassure you of the high quality of our products, and give you added peace of mind should anything happen with these important products
- Choice of flake or nugget models - To suit individual requirements
- Ice is compressed during production - Gives longer lasting cooling for transportation of blood and organs
- Colour: grey
- Climate Class: +5 to +40°C
- Refrigerant Type: R134a
- UK 3 Pin Plug
- Voltage Supply: 230V
- Energy Consumption kw/24h: 7.2
- Freezing Capacity KG per 24H: 65

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Hoshizaki Freestanding Flake/Nugget Ice Machine (85kg/24h)



Product Code: FPD-10233

Dimensions: 800mm H x 640mm W x 600mm D

- Maintenance contract available - Consists of two visits per year at six month intervals to clean, sanitise and check operation, to help prevent breakdown and cut down on staff cleaning times
- 2 years parts and labour warranty - To reassure you of the high quality of our products, and give you added peace of mind should anything happen with these important products
- Choice of flake or nugget models - To suit individual requirements
- Ice is compressed during production - Gives longer lasting cooling for transportation of blood and organs
- Colour: grey
- Climate Class: +5 to +40°C
- Refrigerant Type: R134a
- UK 3 Pin Plug
- Voltage Supply: 230V
- Energy Consumption kw/24h: 7.2
- Freezing Capacity KG per 24H: 85

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Hoshizaki Freestanding Flake/Nugget Ice Machine (125kg/24h)



Product Code: FPD-10235

Dimensions: 800mm H x 640mm W x 600mm D

- Maintenance contract available - Consists of two visits per year at six month intervals to clean, sanitise and check operation, to help prevent breakdown and cut down on staff cleaning times
- 2 years parts and labour warranty - To reassure you of the high quality of our products, and give you added peace of mind should anything happen with these important products
- Choice of flake or nugget models - To suit individual requirements
- Ice is compressed during production - Gives longer lasting cooling for transportation of blood and organs
- Colour: grey
- Climate Class: +5 to +40°C
- Refrigerant Type: R134a
- UK 3 Pin Plug
- Voltage Supply: 230V
- Energy Consumption kw/24h: 12.48
- Freezing Capacity KG per 24H: 125

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Hoshizaki Freestanding Flake/Nugget Ice Machine (125kg/24h)



Product Code: FPD-10237

Dimensions: 800mm H x 940mm W x 600mm D

- Maintenance contract available - Consists of two visits per year at six month intervals to clean, sanitise and check operation, to help prevent breakdown and cut down on staff cleaning times
- 2 years parts and labour warranty - To reassure you of the high quality of our products, and give you added peace of mind should anything happen with these important products
- Choice of flake or nugget models - To suit individual requirements
- Ice is compressed during production - Gives longer lasting cooling for transportation of blood and organs
- Colour: grey
- Climate Class: +5 to +40°C
- Refrigerant Type: R134a
- UK 3 Pin Plug
- Voltage Supply: 230V
- Energy Consumption kw/24h: 12.48
- Freezing Capacity KG per 24H: 125

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Hoshizaki Freestanding Flake/Nugget Ice Machine (110kg/24h)



Product Code: FPD-10238

Dimensions: 800mm H x 940mm W x 600mm D

- Maintenance contract available - Consists of two visits per year at six month intervals to clean, sanitise and check operation, to help prevent breakdown and cut down on staff cleaning times
- 2 years parts and labour warranty - To reassure you of the high quality of our products, and give you added peace of mind should anything happen with these important products
- Choice of flake or nugget models - To suit individual requirements
- Ice is compressed during production - Gives longer lasting cooling for transportation of blood and organs
- Colour: grey
- Climate Class: +5 to +40°C
- Refrigerant Type: R134a
- UK 3 Pin Plug
- Voltage Supply: 230V
- Energy Consumption kw/24h: 12.48
- Freezing Capacity KG per 24H: 110

Delivery

Lead Time: Please check production lead times with our sales.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

How To Order

NHS Purchase Order, MOD Purchase Order, Local Authority Purchase Order

All official government, NHS, MOD, Police, Educational and Local Authority purchase orders are accepted. Fax your purchase order to 0115 727 0123 or email to info@pharmacy-fridge.co.uk.

Bank Transfer

Sent bank transfers to:

Payee	FPD Group Limited		
Bank	HSBC Bank plc		
Sortcode	401903	IBAN	GB97MIDL40190302074613
Account Number	02074613	Swift	MIDLGB2104A

If possible send confirmation of your transaction, to speed up the order process.

Cheques

Make cheques payable to FPD Group Limited and post to:

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Credit/Debit Card

Payment from UK based businesses can be taken over the telephone or processed via the website [shopping cart](#). Customers with issued quotations may process their order online using our [Purchase Code](#) web page.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Delivery

48 Hour delivery is available on most refrigerators through-out the UK mainland, larger units may require a longer lead time, please contact our sales staff whilst ordering.

Approximate Delivery lead times can be confirmed during our quotation process, typical delivery time is 10-14 days from receipt of payment.

**Drivers will not carry goods into buildings.

Goods normally arrive by truck on a pallet and are delivered Kerbside. Please note you may need lifting equipment or a number of strong workers to handle the order when it arrives. If you have any special requirements please [contact us](#).

Report any damages or shortages immediately, please read delivery notes carefully and use contact numbers provided or call FPD Group Ltd.

For more detail refer to [FPD Group Ltd Terms and Conditions](#) documentation.

**Large Products will require a two person delivery and installation, this is built into the price, the larger fridges also may require a longer delivery time of 3-4 weeks please check whilst ordering.

48 HOUR
DELIVERY
AVAILABLE

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

FPDGROUP[™]

Pharmacy Fridge Exports

Since 1961 FPD Group Ltd's family have been supplying medical equipment. Our extensive experience and technical knowledge now concentrates on providing clients with a quality selection of products for hospitals and healthcare organisations worldwide. We receive export enquiries and orders from all parts of the world and have established customers in UK, Europe, Africa and Asia. Export business activities are supervised by either our UK or Asia offices.

Our client base ranges from overseas ministry of health, defence, hospitals, surgeries, clinics and doctors to aid agencies and disaster relief. FPD Group Ltd export almost any medical product listed on our websites.

FPD estimators have the expertise to supply a single item or to fully equip many types of healthcare facilities. As well as the initial supply of equipment we also offer room design services using our in-house Autodesk Revit Architecture 2013 CAD systems.

** Export items are normally subject to minimum order values, additional protective packing charges and export shipping costs. Prices shown are Ex Works unless otherwise agreed in writing.*



FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

About Pharmacy Fridge

Pharmacy Fridge supply a wide range of fridges meeting the specifications of the healthcare industry. To buy simply add your required product to the shopping cart and submit your order.

For dispensary, surgery, ward and general use for the storage of drugs and vaccines under controlled conditions in line with the "Code of Ethics and Standards" of the Royal Pharmaceutical Society of Great Britain. Where necessary to achieve the required temperature control internal air-circulating fans are fitted.

Temperature display, alarms, max/min temperature recording, and door locks are standard and all models are CFC, HCFC and Ammonia free.

Are you storing your vaccines correctly? Recent research has shown an alarming level of use of domestic refrigerators to store vaccines. Did you know that not only can this damage the effectiveness of the vaccines it also poses a security risk for medicines that should be locked away?

FPD Group Ltd

FPD Group Limited is a family business with over 50 years experience in the supply and installation of a wide range of Healthcare Equipment, Medical Refrigeration and Hygienic Storage Furniture. Our customer base includes organisations such as NHS Trusts, PCT's, HMP Prison service, Ministry of Defence, GP Practices, Pharmacy Shops and Dispensaries, Universities, Schools, Dental Practices, Laboratories, Veterinary Practices, Care Homes and Private Buyers. We also export to the Irish Republic, Asia, United Arab Emirates and Qatar. In the UK the Company is based in Exeter, Devon with an additional office in Asia.

Our product range includes pharmacy storage solutions, refrigeration, controlled drugs storage, medical trolleys, clinical waste bins, medical lighting solutions, medical and healthcare furniture ranges and a number of bespoke solutions on request.

FPD Group has an in-house CAD visualisation department and accepts most file formats from your architects' plans, please ask for a quotation.

We are fully committed to the provision of a comprehensive range of professional products, services and furniture. Our philosophy is to provide a complete and appropriate solution to our customers' requirements.

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk

Contact Pharmacy Fridge

FPD Group Ltd have over 40 years experience in the healthcare and shopfitting industry. We offer consultancy, design and installation services. Along with our Pharmacy Shelving and Storage solutions we also provide a wide range of medical equipment and healthcare furniture. Please visit www.fpd-group.co.uk for more information on our medical products.

For general enquiries complete the form below and we will contact you ASAP, please provide as much information as possible including dimensions where required. Please email any plans or photos of installations directly to sales@medicine-fridge.co.uk.

Name	<input type="text"/>	Address	<input type="text"/>
Company	<input type="text"/>	...	<input type="text"/>
Telephone	<input type="text"/>	City	<input type="text"/>
Fax	<input type="text"/>	State	<input type="text"/>
Email	<input type="text"/>	PostalCode	<input type="text"/>
		Country	<input type="text"/>
Comments	<input type="text"/>		

Submit Enquiry

Sales

Tel: 0113 350 8696
Fax: 0115 727 0123
Email: sales@pharmacy-shelving.co.uk

Payments & Billing

Tel: 01392 927 567

Address

FPD Group Limited
The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

FPD Group Limited

The Innovation Centre
University of Exeter, Rennes Drive
Exeter, EX4 4RN

Tel 0113 350 8696
Fax 0115 727 0123
Email info@pharmacy-fridge.co.uk